



2023-2024

Catalogue







# Preface

Dear Customer

We proudly present our new catalogue to you. That is right, MVS In Motion is in a new jacket! Together, we will unwrap our new branding for the coming years on Medica Trade Fair 2023. This move follows the further growth to professionalism and will help us to enlarge our brand awareness in the world of physiotherapy and rehabilitation.

The birth of this new catalogue goes hand in hand with the move to our new offices in Roeselare. Our new offices are the closure of a year filled with construction and tight deadlines. We would also like to invite you, as our customers, to stop by and visit us at our new location. We are delighted to show you our new work environment, our professional showroom and dynamic team.

In order to keep up pace with the fast global changes, we once again have some new products to present to you. Check out our novelties section on page 14 to discover them all.

Of course, all of this would not be possible without you and without our fantastic team. A special thanks to everyone associated with MVS In Motion.

The next step will be the transformation of our website, which will have the same branding as the catalogue. Just bear with us a moment longer, we will inform you about the launch soon.

We kindly invite you to pay us a visit at our booth at Medica 2023 and Arab Health 2024 to celebrate our complete transformation. Let's flip the page together to a new era of MVS In Motion.

Together, we will stretch our limits, fiercely and with a lot of fun!

Sincerely yours,

Evelyn and Tom





# Table of contents

2	Preface
3	Table of contents
7	Brands
9	MVS In Motion
11	Moves
13	Mambo Max
14	New products
<b>19</b>	<b>O1 Resistive</b>
22	Moves band
23	Moves latex-free band
24	Moves tube
25	Moves wide loop
	Moves loop
26	Moves fabric elastic loops
28	Moves FIT band
28	Moves FIT tube
29	Moves FIT loop
	Moves FIT superloop
30	Moves bar
	Moves ring tube
31	Multi-Grip™ Exerciser
	Moves ankleciser
32	Moves shoulder rope pulley
	Door anchor   notch handles
	Door anchor   wooden handles
	Door anchor   hand support
	Metal bracket
33	Value pack
	Moves shoulder tube pulley
	Moves knee rope pulley
34	Moves accessories
<b>35</b>	<b>O2 Hand</b>
36	Moves Flex-Ion
37	Digi-Flex®
	Digi-Extend®
38	Handmaster Plus™
	Prohands® series
39	VariGrip Therapy
	VariGrip Sport
	Constant Force X-Tend
40	Power-Web® Combo
	Power-Web® Senior
	Power-Web® Junior
41	Power-Web® Flex-Grip
	Adjustable hand grip
	Foam hand grip
42	Moves squeeze balls & eggs



43	Moves wrist exerciser
44	Moves comfort putty
45	Theraflex™ putty
	Therapie Sand
<b>46</b>	<b>O3 Upper &amp; lower body</b>
47	Moves oxycycle 1 - active magnetic
	Moves oxycycle 2 - passive
48	Moves oxycycle 3 - active & passive
	Moves oxycycle pedals
49	CuraMotion
50	Moves foldable pedal exerciser
	Moves standard pedal exerciser   premium chrome
	Moves standard pedal exerciser K/D
51	Fit stretch
	Moves Wooden incline board
	Moves wooden incline board - fix
52	Mambo Max adjustable wooden incline board
	Plastic incline board
<b>53</b>	<b>O4 Wellness</b>
54	Moves mini massage gun
55	Triggerpoint rollers
57	Moves IASTM tools
58	Massage tools
59	Massage ball, roll & ring
	Sport tape
	Kinesiology tape
60	Sitting cushions
62	Adjustable ball chair
	Pilates
	Foam Roller
	Soft-Over-Ball
	Ring
64	Stepping stones
	Yoga
<b>65</b>	<b>O5 Fitness</b>
66	Gym mat
67	Aerobic step
	Fabric elastic band
68	Balance board
	Wooden balance board
69	Balance board freeman
	Body twister
	Ekilibre
70	Balance trainer
	Balance pad
71	Balance pad pro
	AB gym ball
72	AB peanut ball
	Dynadome
73	Ball accessories



74

75

76

77

## O6 Strength

SoftMeds

Coated dumbbells

Medicine ball

Kettlebell

Wrist & ankle weights

Premium wrist & ankle weights

Thumb lock wrist weight

78

79

80

81

82

## O7 Thermal

Para Pro - paraffin bath

Coming soon: Hydrocollator and moist heat packs

Hot/cold pack

PearlPack

PearlPack knee/elbow

PearlPack eye mask

PearlPack back wrap

PearlPack sports pack + strap

83

84

Roll-on pack

Instant Heat Pack

85

86

## O8 Evaluation

MicroFET®2 muscle test dynamometer

Fingertip pulse oximeter

Junior oximeter

Senior oximeter

87

Blood pressure monitor

Infrared thermometer

Ear & forehead thermometer

Shotgun thermometer

88

Pedi-scolimeter

Neurological reflex instruments

89

90

## Saehan

Hydraulic hand dynamometer

Hydraulic 3-piece hand evaluation kit

Hydraulic 7-piece hand evaluation kit

91

Digital hand dynamometer with data transfer software - DHD-1

Digital hand dynamometer with G-STAR™ software - DHD-3

Squeeze dynamometer

92

Smedley hand dynamometer

Back-leg-chest dynamometer

Pinch gauge

93

Plastic goniometer

Stainless steel goniometer

94

Finger goniometer

Goniometer kit

Tuning fork

95

Fat caliper

96

97

## Baseline®

BIMS™ digital (Load-Cell) 5-position dynamometers

Functional model



	Clinic model
	Deluxe model
99	5-position hydraulic hand grip dynamometers
	Pneumatic squeeze dynamometers
100	Smedley spring dynamometer
	Electronic Smedley dynamometer
	Mechanical pinch gauge
101	Hydraulic pinch gauge
	5-position hydraulic pinch gauge
	Electronic pinch gauge
102	7-piece hand sets
103	3-piece hand sets
104	Push-pull dynamometer
	Back-leg-chest dynamometer
	Finger (digit) goniometer
105	Plastic finger goniometer
	180° 20 cm goniometer
	180° 35 cm goniometer
	360° 35 cm goniometer
106	180° finger (small joint) goniometer
	6-piece goniometer set
	XTender® extendable goniometer
107	Lafayette goniometer
	Plastic goniometers
108	Digital goniometer
	Touch-Test® 2-point discriminator
	Discrim-A-Gon® 2-point discriminator
	2-point discriminators
109	Sit 'n Reach® box test
	Finger flexion gauge
	Gulick measurement tape
	Finger circumference gauge
	Circumference tape
110	Clinical monofilament sensory evaluators
111	Tuning Forks
	Pinwheel
112	Bubble® inclinometer
	AcuAngle® inclinometers
	Digital inclinometer
	Scoliometer
113	Economy skinfold caliper
	Skinfold calipers
	BMI body fat scale
114	Private label
116	After-sales
118	Terms & conditions



# Brands

## Our own brands



---

## Master distributor for









MVS In Motion is a dynamic and quality-driven B2B company, with people dedicated to excellence in everything they do. We're a master distributor, specializing in physiotherapy, rehabilitation, sports, and fitness products. In 2021 the company was acquired by the new owners, Evelyn Deceuninck and Tom Vanhecke.

In 2023 we relocated to Roeselare (Belgium). True to our character, we completely reinvented ourselves. Now, we can proudly present the new and improved MVS In Motion.



#### ... **EXPERTISE**

MVS In Motion is that one-stop shop that helps you to grow and succeed. Apart from offering valuable products, we support our customers with unique knowledge and expertise. Our basis is a good structure and well-oiled organization that efficiently responds to every question.



#### ... **PASSION**

We make a difference through our passion. MVS In Motion is committed and does not give up. We are continuously taking steps to deliver the best.



#### ... **DARING**

MVS In Motion is constantly looking out for interesting and relevant opportunities. We do this to keep up with the market needs.



#### ... **VIBRANT**

Innovating, evolving, and excelling is in our DNA. We are constantly on the move and we bring out the best in ourselves as well as in our customers. We are vibrant and encouraging, thanks to the solid foundation and the strong MVS team. You can always rely on us.



#### ... **INSIGHT**

Insight into your needs and wishes is crucial. That is why we keep improving our wide range of products. We are energetic, powerful, and we like to tackle anything head-on.



**moves** 



**STRETCHING LIMITS**



Moves is one of the own brands of MVS In Motion. It is a top-notch brand for people who want honest and terrific products. Moves helps you to become a better version of yourself. It helps you to stretch your limits.

Under the Moves brand, we also offer Moves FIT products. These products have a unique colour sequence which distinguishes us from the known colours available on the market.



## ... **ESSENCE**

Moves is there for you. Our focus and high standards is what makes us able to assist you with essential products and specific expertise.



## ... **HONEST**

We are ready for each and every request for help. We deliberately opt for a wide but specialised product range. These products have long proven their worth. Honest products that excel in quality and convenience.



## ... **EXPERTISE**

Moves stands for expertise and professional know-how. A quality-driven brand that guarantees optimal support. Moves products help during your rehabilitation or movement optimization. They help you to achieve your goals and succeed.



## ... **EVOLUTION**

With the aid of Moves, you move forward, step by step, towards your desired progression.

There are products for every phase, for every treatment. Moves is there to assist your evolution. You can bring your patient to the next level. With Moves, we help you to reach the top.



## ... **MEANINGFUL**

Moves keeps reinventing and improving its products. This brand has the ambition to do better each time. Moves makes the difference. We offer the best selection of products, that supports and eases the commitment and discipline of your patients.



mambo  
max

Fit.  
Fun.  
Fierce.





Mambo Max is the other brand of MVS In Motion. It's a light-hearted and fun brand that assists anyone in exercising with joy. Mambo Max lights the fire to become a healthier and happier version of yourself. Get fit, become fierce, but most importantly: have a lot of fun!



... **ENERGY**

Mambo Max brings balance into your life. This brand inspires you to take action, to work out, to move. It even assists you in relaxation and massage. It is the perfect partner for a healthy life.



... **FIT**

A body that is in shape, is a healthy body. Mambo Max offers a complete assortment to keep your body fit. It helps you to build power, resilience, flexibility, and balance. Together with Mambo Max, you move forward, with the right balance.



... **FUN**

Mambo Max is a symbol for fun. Our colourful range of products emphasises this. The goal is to give you that optimal, fun workout. Exercising becomes a true pleasure. You have to admit, working on yourself doesn't have to feel like working, does it?



... **ATTRACTIVE**

Mambo Max is here for everyone. Young and old, skilled or newbie, for people that want to see results. We make sports and exercising attractive. The happy and appealing colours guarantee a light-hearted, playful range.



# New products



## Fabric elastic loops set

This trendy loop set is the perfect solution for building your muscle power in an easy and effective way. Due to the anti-slip pattern, they stay exactly where they need to be and you don't have to worry about them slipping off.

- Light grey: light resistance
- Dark grey: medium resistance
- Black: heavy resistance

### SPECIFICATIONS

- ✓ Comes in a mesh bag
- 👁️ Set contains 3 different resistances: light | medium | heavy

Description	Reference code	Unit
Fabric Elastic Loops	01-704001	set of 3



## Kinesiology tape 35m roll MD CE

Moves Kinesiology Tape is used **to facilitate motion, inhibit pain and improve circulation** while still supporting the muscles and joints. Until recently, this tape was only available in 5m rolls. We have heard your requests carefully and now the Moves Kinesiology tape in tan is also available in a 35m roll.

### SPECIFICATIONS

- 📏 Dimensions:
  - 5 cm x 35m



Colour	Description	Dimensions	Reference code	Unit
🟡	Moves Kinesiology Tape   Tan	5 cm x 35 m	04-060105	1 pc




## MicroFET™2 muscle test dynamometer



The wireless microFET™2 Digital Handheld Dynamometer is an accurate, portable Force Evaluation and Testing (FET) device. It is designed for taking objective, reliable and quantifiable **muscle testing measurements**. The microFET™2 aids in diagnosis, prognosis and treatment of neuromuscular disorders.

- Measurement range 0-136 kg (0-300 lb) force
- Easy-to-read LCD Displays shows peak force and elapsed time
- Selectable units of measure: pounds (lbs.), Newtons (N), or kilogram-force (kgf)
- Accuracy within 1% of the reading
- Two threshold settings for muscle testing: Low 0,4 to 136 kg (0,8 to 300 lb) in 0,05 kg (0,1 lb) increments and High threshold 1,4 to 136 kg (3,0 to 300 lb) in 0,05 kg (0,1 lb) increments
- Stores up to 30 tests
- Uses rechargeable lithium ion rechargeable battery
- Self-activating "sleep" mode after 3 minutes to extend battery life

Description	Reference code	Unit
MicroFET™2 Muscle Test Dynamometer	12-O381W	1 pc



## Premium wrist & ankle weights

Complete your gym set and enhance your work out with these trendy wrist and ankle weights. With the premium Mambo Max wrist and ankle weights you can easily train your arm and leg muscles. You can build up the weights or simply challenge yourself during exercising.

Available in 3 different weights, each with a particular colour:

- Light grey: 0,5 kg
- Grey: 1 kg
- Anthracite: 1,5 kg

### SPECIFICATIONS

- Available in 3 different weights:  
0,5 kg | 1 kg | 1,5 kg
- Pair of the chosen weight comes in a colourbox

Description	Reference code	Unit
Wrist & ankle weights   0,5 kg   Light grey	06-O60101	pair
Wrist & ankle weights   1 kg   Grey	06-O60102	pair
Wrist & ankle weights   1,5 kg   Anthracite	06-O60103	pair





## Adjustable ball chair

The Mambo Max adjustable ball chair offers you the ideal way to upgrade your work setup. The ergonomic back support assists with spine alignment, it strengthens the core, and it corrects your sitting posture.

- Instability-induced motion helps to firm your core, improve circulation and posture
- Chair height ranges from 51 cm to 58 cm (20" to 23") to fit most standard desks
- Stable, 360° swivel base allows for easy movement
- Lightweight and compact
- Lockable easy glide caster wheels
- Anti-burst, removable oval-shaped ball seat, different from other ball chairs
- Max weight: 127 kg (280 lbs)



Colour	Description	Reference code	Unit
	Mambo Max Adjustable Ball Chair Black/Grey	O4-O2O4O1	1 pc


## Xtra comfort gym mat

Mambo Max Gym Mats provide a soft touch and comfortable exercise environment. They are perfect for gym & fitness. They have a non-slip, well-cushioned surface with a unique closed cell structure. The mats are equipped with plastic rings that are 40 cm apart so they can be suspended while not in use. Mambo Max also offers a Gym Mat Wall Hanger that can hold up to 10 mats. They are **latex-free** and can be used in water. Each mat comes with a handy **storage/carry strap**.

The Mambo Max xtra comfort gym mat is now also **available in black**.

### SPECIFICATIONS

 measurements: 180 x 60 x 1,5 cm



Colour	Reference code	Description	Unit
	O5-O3O1O4	Mambo Max Xtra Comfort Gym Mat   180 x 60 x 1,5 cm   Rings   Strap	1 pc



Multi-Filament Pen

MD

CE

The multi-filament pen offers a non-invasive evaluation of the diminished protective sensation. The pen contains 4 nylon monofilaments, each representing a unique amount of force which is marked by colour and size. Simply slide down the lever of the desired monofilament and use as intended. This pen provides a useful set of monofilaments in a fun, compact, and protected form.

Reference code	Description	Unit
12-1370	Baseline Multi-Filament Pen	1 pc



Therapie Sand

MD

CE

Therapie Sand is an innovative material that provides a soothing tactile-kinesthetic experience and is equally popular among both children and adults. The sand promotes fine motor skills and sensitivity, stimulates circulation to relieve everyday stress, and restores mobility and overall strengthening of the wrists and ankles.

It is primarily used in therapy to promote sensitivity and optimize mobility. It improves overall strength and hand-foot coordination. Additionally, the Therapie Sand is an excellent educational aid, particularly in the fields of occupational therapy and physiotherapy, for the purpose of stimulating creativity and development.

SPECIFICATIONS

- Q

Material: Latex-free sand
- 🛒

Available in 2 different weights: 454 g - 2 kg



Colour	Reference code	Description	Unit
<div>●</div>	O2-24O4O1	Therapie Sand   454 g   Grey	1 pc
	O2-24O5O1	Therapie Sand   2 kg   Grey	1 pc





### Adjustable wooden incline board

The Mambo Max adjustable wooden incline board completes your workout setup. It is crafted from durable poplar wood and has a multi-layered board. This board can be adjusted in 4 different angles and conveniently folds for easy storage. Its slip-resistant stepping surface ensures stability during workouts. Experience enhanced stretching, rehabilitation, and strength training with this versatile incline board, perfect for targeting various muscle groups and improving flexibility and balance.

#### SPECIFICATIONS

- ↔ Measurements: 46 x 38 x 5,8 cm
- Q Material: poplar wood
- ∠ Inclination: adjustable to 10°, 20°, 30°, or 40° angles



Reference code	Description	Unit
O3-O40002	Mambo Max adjustable wooden incline board	1 pc





# O1 Resistive





## One-stop-shop for all your bands and tubes

MVS In Motion is your one-stop-shop for all resistive exercise products. With our premium quality brand Moves, we offer both original colour and FIT colour sequence products. In addition, we offer a wide range of Latex-Free resistive exercise bands in original colours. Our tubes and loops are also available in the original and FIT colour sequence.

The latex bands, tubes and (super)loops are made of high quality latex from huge rubber tree plantations in Malaysia. All products have an extraordinary elasticity and linearity. They are the ideal partner for your physiotherapy trainings or intensive work-outs, because:

- ✓ our bands have the patented Snap-Stop® feature
- ✓ our tubes go through a dipped production process
- ✓ our (super)loops are equipped with a multi-layered structure

The different resistances allow you to use Moves and Moves FIT products for rehabilitation as well as for fitness purposes.

Moves Latex-Free is made from TPE, with the same high quality standards as the regular Moves bands. Never before were the characteristics from Latex-Free bands so close to those from their latex counterpart.

- ✓ **Moves Band:** 8 resistances (Tan - Gold) & 5 lengths (1,5m - 2,5m - 5,5m - 22,5m - 45,5m)
- ✓ **Moves FIT Band:** 6 resistances (Peach - Grey) & 4 lengths (1,5m - 2,5m - 5,5m - 45,5m)
- ✓ **Moves Latex-Free Band:** 5 resistances (Yellow - Black) & 3 lengths (1,5m - 5,5m - 45,5m)

Original colours		Same resistance	FIT colours	
	Tan   Extra light			Peach   Light
	Yellow   Light			Orange   Medium
	Red   Medium			Lime green   Heavy
	Green   Heavy			Blueberry   Extra heavy
	Blue   Extra heavy			Plum   Special heavy
	Black   Special heavy			Grey   Super Heavy
	Silver   Super heavy			
	Gold   Ultra heavy			

- Well-known colours used by physiotherapists all over the world
- Moves brings premium quality at an affordable price range

- FIT colours are more trendy and youthful than the original colours
- Sell Moves FIT next to Moves to widen customer base and generate more sales





## What is Snap-Stop®?

Moves and Moves FIT Bands are unique in their field. Thanks to our extensive research and development, we have achieved an innovative production process: **the Snap-Stop® technology**.

This distinctive feature will avoid sudden snapping when the band is exposed to high forces, even if there is a cut or puncture! This makes exercising with our bands much safer.



SNAP-STOP® FOR SAFER RESISTIVE EXERCISING



## Elongation

Moves Bands can be stretched up to 300% elongation and Moves Loop Bands even up to 400% elongation! Therefore they are your perfect companion for a rehabilitation programme or a full body work-out. Most competitor's bands are not to be stretched more than 3 times their resting length, while a Moves Band easily achieves more, maintaining its amazing elasticity.

The high-quality material enables you to have a very wide array of resistances and possibilities.

### Note:

50% elongation = 1,5 x lay-flat length

100% elongation = 2 x lay-flat length

150% elongation = 2,5 x lay-flat length

300% elongation = 4 x lay-flat length

400% elongation = 5 x lay-flat length



## Band MD CE

Premium quality, powder-free, latex resistive exercise bands. Tear-resistive thanks to the Snap-Stop® fabrication process. The Moves Band provides positive as well as negative forces on the muscles. Furthermore it improves the range of motion, the strength and the cooperation of the muscle groups.









### SPECIFICATIONS










 Lay-flat-length: available in 5 different lengths

 Width: 127mm

 Resistance: available in 8 different levels



Total Kg / Elongation %						
	50%	100%	150%	200%	250%	300%
	0,8	1,1	1,3	1,6	1,9	2,2
	0,9	1,4	1,7	2,0	2,3	2,7
	1,3	1,8	2,2	2,6	3,1	3,7
	1,6	2,3	2,8	3,3	3,9	4,6
	1,9	2,7	3,3	4,0	4,7	5,5
	2,5	3,6	4,5	5,3	6,2	7,3
	3,5	5,0	6,1	7,3	8,6	10,1
	3,8	5,4	6,7	7,9	9,4	11,0

Colour	Resistance	1.5m Dispenser 25 pcs	1.5m 10-pack	2.5m 10-pack	5.5m Box Small Roll	22.5m Box Medium Roll	45.5m Large Roll
	Tan   Extra Light	n/a	n/a	n/a	01-100501	n/a	01-104501
	Yellow   Light	01-100102	01-101502	01-100202	01-100502	01-102202	01-104502
	Red   Medium	01-100103	01-101503	01-100203	01-100503	01-102203	01-104503
	Green   Heavy	01-100104	01-101504	01-100204	01-100504	01-102204	01-104504
	Blue   Extra Heavy	01-100105	01-101505	01-100205	01-100505	01-102205	01-104505
	Black   Special Heavy	01-100106	01-101506	01-100206	01-100506	01-102206	01-104506
	Silver   Super Heavy	n/a	01-101507	01-100207	01-100507	01-102207	n/a
	Gold   Ultra Heavy	n/a	n/a	n/a	01-100508	01-102208	n/a
	Set of 5 resistances	n/a	01-101511 (5-pack)	01-100211 (5-pack)	n/a	n/a	n/a

### PACKAGING

- ✓ Dispenser 1,5m
- ✓ Polybag: 10 individually packed pcs of the same resistance and length
- ✓ Box with one roll: available in 3 different lengths (5,5m / 22,5m / 45,5m)



## Latex-free band



Latex-Free exercise band with characteristics that are extremely close to its latex counterpart. Exceptional elasticity and 'resistance vs. elongation' linearity.

Moves Latex-Free Bands start from superior thermoplastic elastomer (TPE), so it's non-toxic and completely PVC and phthalate free.

### SPECIFICATIONS






 Lay-flat-length: available in 3 different lengths

 Width: 100mm

 Resistance: available in 5 different levels






#### Total Kg / Elongation %

	50%	100%	150%	200%	250%	300%
	1,9	2,4	2,7	3,2	3,8	4,6
	2,2	2,8	3,2	3,7	4,5	5,3
	2,5	3,1	3,7	4,3	5,1	6,1
	3,2	4,1	4,8	5,6	6,6	7,9
	4,0	5,1	6,0	7,0	8,3	9,9

Colour	Resistance	1.5m 10-pack	5,5m Small roll	45,5m Large roll
	Yellow   Light	O1-1215O2	O1-1205O2	O1-1245O2
	Red   Medium	O1-1215O3	O1-1205O3	O1-1245O3
	Green   Heavy	O1-1215O4	O1-1205O4	O1-1245O4
	Blue   Extra Heavy	O1-1215O5	O1-1205O5	O1-1245O5
	Black   Special Heavy	O1-1215O6	O1-1205O6	O1-1245O6
    	Set of 5 resistances	O1-121511 (5-pack)	n/a	n/a

### PACKAGING

-  Polybag: 10 individually packed pcs of the same resistance and length
-  Box with one small roll: 5,5m
-  Box with one large roll: 45,5m



**Tube** MD CE

Our high-standard Moves Tubes are created by an intensive dipping process. The multiple layers guarantee a strong and safe training product. Extruded tubes are more susceptible to damage and can easily snap. So, our multi-layered dipped tubes are a superior and durable product. Train all the muscles in your body and improve your body flexibility.

Be assured of an exceptional strength to elongation linearity from our tubes.



	Total Kg / Elongation %					
	50%	100%	150%	200%	250%	300%
	1,0	1,5	1,9	2,3	2,9	3,6
	1,4	2,0	2,5	2,8	3,4	4,0
	2,4	3,4	4,0	4,9	5,5	6,5
	2,9	3,9	4,7	5,6	5,9	7,4
	3,4	4,4	5,4	6,2	7,1	8,1
	4,0	5,4	6,5	7,7	8,1	9,5

Colour	Resistance	120cm with handles	7,5m roll	30m roll
	Yellow   Light	01-200102	01-200702	01-203002
	Red   Medium	01-200103	01-200703	01-203003
	Green   Heavy	01-200104	01-200704	01-203004
	Blue   Extra Heavy	01-200105	01-200705	01-203005
	Black   Special Heavy	01-200106	01-200706	01-203006
	Silver   Super Heavy	01-200107	01-200707	01-203007

**PACKAGING**

- ✓ Polybag: 120cm tube with soft handles
- ✓ Box with 1 roll of 7,5m
- ✓ Box with roll of 30m | contains maximal 3 separate pcs - minimal length of 1 pc is 4m



COMBINE A ROLL WITH OUR SAFETY HANDLES (P.19) FOR SUPERIOR GRIP







## Wide loop



We offer both regular and wide loops. Wide loops are ideal for spreading forces more evenly on your body parts. They are high quality, cost effective and easy to use. The Moves Wide Loops are multi-layered and available in different resistances, thus satisfying all users.

### SPECIFICATIONS

Lay-flat-length: 30cm

Width: 7,5cm

Colour	Resistance	Reference code 10-pack
	Tan   Extra light	01-100071
	Yellow   Light	01-100072
	Red   Medium	01-100073
	Green   Heavy	01-100074
	Blue   Extra heavy	01-100075
	Black   Special heavy	01-100076



### PACKAGING

1 unit: 10 individually packed pcs of the same colour/resistance



## Loop



Moves Loops are extremely portable and unbelievably versatile products. They are very effective instruments for stretching and strengthening all muscle groups. The Moves Loops are multi-layered, which makes them thicker and more durable. The benefit of these thicker loops is that they won't roll up during exercising, which avoids the uncomfortable feeling from thinner alternatives.

### SPECIFICATIONS

Lay-flat-length: 30cm

Width: 2,5cm

Contains natural latex rubber

### PACKAGING

Polybag: 10 individual pcs of the same colour/resistance.

Mesh bag: Set of 4 pcs of the most popular resistances

Colour	Resistance	Reference code 10-pack
	Tan   Extra light	01-100001
	Yellow   Light	01-100002
	Red   Medium	01-100003
	Green   Heavy	01-100004
	Blue   Extra heavy	01-100005
	Black   Special heavy	01-100006
	Set of 4 resistances	01-100011



Total Kg / Elongation %						
	50%	100%	150%	200%	250%	300%
	1,9	2,8	3,4	4,0	4,7	5,3
	2,4	3,5	4,3	5,0	5,8	6,7
	3,1	4,5	5,6	6,6	7,6	8,7
	3,8	5,6	6,9	8,1	9,3	10,7
	4,8	7,0	8,6	10,1	11,7	13,3
	6,0	8,7	10,7	12,6	14,6	16,7



**New: Fabric elastic loops set**

This trendy loop set is the perfect solution for building your muscle power in an easy and effective way. Due to the anti-slip pattern, they stay exactly where they need to be and you don't have to worry about them slipping off.

- Light grey: light resistance
- Dark grey: medium resistance
- Black: heavy resistance

**SPECIFICATIONS**

- ✓ Comes in a mesh bag
- 👁 Set contains 3 different resistances:  
light | medium | heavy

Description	Reference code	Unit
Fabric Elastic Loops	OI-7O4OOI	Set of 3







## FIT band MD CE



Moves FIT Bands are powder-free, latex resistive exercise bands with the same premium quality as the original Moves Bands. The difference is purely the alternative, more youthful colour sequence.

### SPECIFICATIONS




-  Length: available in 5 different lengths
-  Resistance: available in 6 different levels



	Total Kg / Elongation %					
	50%	100%	150%	200%	250%	300%
	0,9	1,4	1,7	2,0	2,3	2,7
	1,3	1,8	2,2	2,6	3,1	3,7
	1,6	2,3	2,8	3,3	3,9	4,6
	1,9	2,7	3,3	4,4	4,7	5,5
	2,5	3,6	4,5	5,3	6,2	7,3
	3,5	5,0	6,1	7,3	8,6	10,1

Colour	Resistance	1.5m 10-pack	2.5m 10-pack	5.5m Small roll	22.5m Medium roll	45.5m Large roll
	Peach   Light	01-111502	01-110202	01-110502	n/a	01-114502
	Orange   Medium	01-111503	01-110203	01-110503	n/a	01-114503
	Lime green   Heavy	01-111504	01-110204	01-110504	n/a	01-114504
	Blueberry   Extra Heavy	01-111505	01-110205	01-110505	n/a	01-114505
	Plum   Special Heavy	01-111506	01-110206	01-110506	n/a	01-114506
	Grey   Super Heavy	01-111507	01-110207	01-110507	01-112207	n/a (5-pack)
	Set of 5 resistances	01-111511 (5-pack)	01-110211 (5-pack)	n/a	n/a	n/a

### PACKAGING

-  5-pack: Set of 5 resistances (peach to plum) pre-cut in pieces in one pack
-  10-pack: 10 pre-cut pieces of the same resistance in one pack
-  Box with one roll: available in 5,5m / 22,5m / 45,5m







**FIT tube** MD CE

Just like the Moves Tubes, the Moves FIT Tubes are also produced using an intensive dipping process. Therefore, these multi-layered tubes ensure a safe and strong training. This is a big advantage over extruded tubes, which are more susceptible to damage and can easily snap. Get these hip & trendy FIT colours now!


**SPECIFICATIONS**

 Resistance: available in 6 different levels

Total Kg / Elongation %						
	50%	100%	150%	200%	250%	300%
	1,0	1,5	1,9	2,3	2,9	3,6
	1,4	2,0	2,5	2,8	3,4	4,0
	2,4	3,4	4,0	4,9	5,5	6,5
	2,9	3,9	4,7	5,6	5,9	7,4
	3,4	4,4	5,4	6,2	7,1	8,1
	4,0	5,4	6,5	7,7	8,1	9,5

Colour	Resistance	30m roll
	Peach   Light	OI-213002
	Orange   Medium	OI-213003
	Lime green   Heavy	OI-213004
	Blueberry   Extra Heavy	OI-213005
	Plum   Special Heavy	OI-213006
	Grey   Super Heavy	OI-213007

**PACKAGING**

 Dispenser box with one roll: available in 30m



COMBINE A ROLL WITH OUR SAFETY HANDLES (P.26) FOR SUPERIOR GRIP. MOVES TUBES WITH HANDLES ARE USEFUL FOR PHYSIOTHERAPY AND MUSCLE TRAINING.



## FIT loop



Our excellent and durable loops are also available in Moves FIT colours. The Moves FIT Loops are easy to use and very versatile. You can stretch and strengthen all muscle groups.

These Loops are multi-layered which makes them thicker and more durable.

### SPECIFICATIONS

- Lay-flat-length: 30cm
- Width: 2,5cm
- Resistance: available in 4 levels
- Contains natural latex rubber

### PACKAGING

- Polybag: 10 individual pcs of the same colour/resistance
- Mesh bag: Set of 4 pcs of the most popular resistances



Colour	Resistance	Reference code 10-pack
	Peach   Light	01-110002
	Orange   Medium	01-110003
	Lime green   Heavy	01-110004
	Blueberry   Extra heavy	01-110005
	Set of 4 resistances	01-110011

## FIT superloop



With the Moves FIT Superloop you can easily stretch and strengthen different muscle groups. It can be used for rehabilitation, fitness and strength training.

The Moves FIT Superloop is also multi-layered, which makes it a very durable product. Because of its long lay-flat length and circumference, the FIT Superloop allows you to stretch and strengthen several muscle groups at the same time.



### SPECIFICATIONS

- Lay-flat length: 104cm
- Circumference: 208cm
- Resistance: available in 5 different levels

Colour	Resistance	Reference code 1 pc
	Peach   Light	01-110402
	Orange   Medium	01-110403
	Lime green   Heavy	01-110404
	Blueberry   Extra heavy	01-110405
	Plum   Special heavy	01-110406

		Total Kg / Elongation %					
	Width x Thickness in mm	50%	100%	150%	200%	250%	300%
	15 x 3,0	4,3	6,0	7,2	8,4	9,7	11,1
	22 x 4,0	8,5	11,8	14,1	16,5	19,0	21,7
	30 x 4,4	12,7	17,6	21,2	24,7	28,5	32,5
	40 x 4,4	16,9	23,5	28,3	32,9	38,0	43,4
	53 x 4,4	22,4	31,2	37,4	43,6	50,3	57,5




## Bar



Golfer's or tennis elbow? The Moves Bar is the ideal product to treat and resolve this!

This effective product is developed to focus on hand, arm and shoulder problems. This product allows flexion and tension exercises, adduction and abduction exercises and even soft tissue mobilisation. So, far more possibilities than only twisting and bending!



### SPECIFICATIONS

- Length: 31cm
- Width: 4,5cm
- Material: dried natural rubber

### PACKAGING

- Polybag: containing 1 pc

Colour	Resistance	Kg-Force*	Reference code	Unit
	Extra light	1,8	01-600101	1 pc
	Light	3,0	01-600102	1 pc
	Medium	5,0	01-600103	1 pc
	Heavy	7,0	01-600104	1 pc
	Extra heavy	11,0	01-600105	1 pc
	Special heavy	14,0	01-600106	1 pc

\*FORCE APPLIED WHEN BENDING THE BAR TO A U-SHAPE

## Ring tube



The **Moves Cuff-ring Tube** is a versatile exercise tool to tone and strengthen all the muscles of the lower body. With velcro closure that can be adjusted to all sizes.

The **Moves 8-Ring Tube** offers lots of possibilities for short and long pull exercises thanks to its centre. Used for intensive training of isolated muscles. Efficient for all muscles of the upper- and lower body.


With the **Moves O-Ring Tube** you get an intensive training of each muscular group. Ideal for abdomen, bottom, hip and leg exercises.



MOVES CUFF-RING TUBE



MOVES 8-RING TUBE



MOVES O-RING TUBE

Colour	Resistance	Reference code			Unit
		Cuff-Ring Tube Loop	Figure 8-Ring Loop	O-Ring Tube Loop	
	Yellow   Light	01-200402	01-200302	01-200202	1 pc
	Red   Medium	01-200403	01-200303	01-200203	1 pc
	Green   Heavy	01-200404	01-200304	01-200204	1 pc
	Blue   Extra heavy	01-200405	01-200305	01-200205	1 pc
	Black   Special heavy	01-200406	01-200306	01-200206	1 pc



## Multi-Grip™ Exerciser



Have unlimited exercise options with the Multi-Grip™ Exerciser. Each exerciser contains both large and small loops that are easy to grip. Use small loops as an anchor or for finger exercises. Use the larger loops around the hands and feet to perform high impact upper and lower body exercises. Change resistance easily by moving 1 loop up or down the exerciser or use a different coloured exerciser out of the 8 different resistance levels.

### SPECIFICATIONS

- Length: 30 yard roll (27,4m)
- Width: 4,5cm
- Material: unscented cloth that can be washed and dried. Does not contain latex or a powder coating.

### PACKAGING

- Roll: 30 yards (27,4m)



Total Kg / Elongation %			Colour	Resistance	Reference code	Unit
100%	200%					
	0,9	2,3		Tan   XX Light	01-703001	1 pc
	1,4	3,4		Yellow   X-Light	01-703002	1 pc
	1,8	4,5		Red   Light	01-703003	1 pc
	2,3	5,7		Green   Medium	01-703004	1 pc
	2,7	6,8		Blue   Heavy	01-703005	1 pc
	3,2	7,9		Black   X-Heavy	01-703006	1 pc
	4,5	11,3		Silver   XX-Heavy	01-703007	1 pc
	5,9	14,7		Gold   XXX-Heavy	01-703008	1 pc

## moves

### Ankleciser

The Moves Ankleciser is the ideal product to exercise your ankles. It's as easy in use as it is effective. The feet slip comfortably into the band which consists of two stitched loops. This progressive exercise system enables the user to choose the desired resistance level to perform the exercises.



### SPECIFICATIONS

- Length: 31cm
- Resistance: available in 4 different levels
- Material: latex-free elastic band

### PACKAGING

- Polybag: 5 pcs of the same colour/resistance
- Set of 4 different colours/resistances



Colour	Resistance	Reference code	Unit
	Yellow   Light	01-500102	5 pcs
	Red   Medium	01-500103	5 pcs
	Green   Heavy	01-500104	5 pcs
	Blue   Extra heavy	01-500105	5 pcs
	Set of 4 resistances	01-500111	4 pcs



## Shoulder rope pulley


The Moves Shoulder Rope Pulley is an economical, easy-to-use device that helps increase range of shoulder motion. It can be used for prevention and rehabilitation of shoulder injuries.

### 1. DOOR ANCHOR | NOTCH HANDLES

**MD** **CE**

This product comes with a pulley, door anchor, nylon cord, stopper system and notch handles. The unique notch handles allow the user to easily adjust the length of the rope through the handle while the revolutionary stopper system is a great aid for consistent and progressive exercising.

Reference code	Unit
01-400101	1 pc



DOOR ANCHOR, NYLON CORD AND STOPPER SYSTEM




NOTCH HANDLE

### 2. DOOR ANCHOR | WOODEN HANDLES

**MD** **CE**

This version of the Moves Shoulder rope pulley comes with a pulley, door anchor, nylon cord, stopper system and wooden handles.

Reference code	Unit
01-400105	1 pc



DOOR ANCHOR, NYLON CORD AND STOPPER SYSTEM




WOODEN HANDLE

### 3. DOOR ANCHOR | HAND SUPPORT


**MD** **CE**

This product comes with a pulley, door anchor, nylon cord, stopper system and detachable hand support. This hand support allows patients without gripping power to use the shoulder pulley.

Reference code	Unit
01-400103	1 pc



DOOR ANCHOR, NYLON CORD AND STOPPER SYSTEM



HAND SUPPORT

### 4. METAL BRACKET


**MD** **CE**

This product comes with a pulley, door anchor, nylon cord, stopper system and metal bracket. This allows you to easily mount the shoulder rope pulley in a doorway.

Reference code	Unit
01-400102	1 pc



METAL BRACKET



DOOR ANCHOR, NYLON CORD AND STOPPER SYSTEM




WOODEN HANDLE



## 5. VALUE PACK

The Moves shoulder rope pulley VALUE PACK comes in a polybag, with notch handles but **without stoppers** on the nylon cord.

Reference code	Unit
01-400104	1 pc



DOOR ANCHOR AND NYLON CORD



NOTCH HANDLE

**moves**

### Shoulder tube pulley


MD C E

Just like the Moves Shoulder Rope Pulley, the Moves Shoulder Tube Pulley helps increase range of shoulder motion. It can be used for prevention and rehabilitation of shoulder injuries.




#### SPECIFICATIONS

 Resistance: available in 3 different levels

#### PACKAGING

 Pack: containing 1 tube of the selected colour/ resistance, soft handles, pulley and door anchor



Colour	Resistance	Reference code	Unit
	Yellow   Light	01-400202	1 pc
	Red   Medium	01-400203	1 pc
	Green   Heavy	01-400204	1 pc


**moves**

### Knee rope pulley

MD C E

The Moves Knee Rope Pulley helps increase range of knee motion. This easy-to-use product can be mounted in a doorway to help you recover from or prevent knee injuries.

While sitting on a chair, put your foot into the ankle strap. A rope starts from the ankle strap, goes through a pulley and ends at a wooden handle designed for a 2-handed grip. Pull your foot carefully upwards by moving the handle towards your abdomen which will generate a rotational movement at your knee.



Reference code	Description	Unit
01-400301	Moves Knee Rope Pulley   Door anchor   2-handed Wooded handle   Ankle strap	1 pc



**Accessories**




SAFETY HANDLES



SOFT HANDLES



DOOR ANCHOR



UNIVERSAL ASSIST



MOVES BAND & TUBING SCISSORS

Reference code	Description	Unit
01-300101	Moves Safety Handles	pair
01-300102	Moves Soft Handles	pair
01-300103	Moves Door Anchor	1 pc
01-300104	Moves Universal Assist	1 pc
01-300106	Moves Band & Tubing scissors	1 pc



# O2 Hand





Flex-Ion



The Moves Flex-Ion is the ultimate rehabilitation product after a hand injury. The spring-loaded buttons isolate each digit so stronger fingers cannot compensate for the weaker ones. This results in a better, faster & more controllable rehab process.








The gradual resistance system of the Moves Flex-Ion helps therapists measure progress. Additionally, the ergonomic base adds comfort and stability.

The Moves brand assures a top-notch product!



SPECIFICATIONS

- Resistance: available in 6 different levels
- Material: display in metal

Colour	Resistance: 1 digit / all digits	Reference code	Unit
	Tan   0,35 kg / 1,1 kg	O2-000101	1 pc
	Yellow   0,7 kg / 2,3 kg	O2-000102	1 pc
	Red   1,4 kg / 4,5 kg	O2-000103	1 pc
	Green   2,3 kg / 7,3 kg	O2-000104	1 pc
	Blue   3,2 kg / 10,4 kg	O2-000105	1 pc
	Black   4,1 kg / 14,1 kg	O2-000106	1 pc
 Set of 6 resistances in 1 display		O2-000111	set of 6 pcs
Empty display		O2-000121	1 pc

PACKAGING

- One Flex-Ion of the selected colour/ resistance
- Set of 6 different colours/resistances (with display)



UNIVERSAL ASSIST



SET WITH DISPLAY



# Digi-Flex®



Digi-Flex® develops isolated finger strength, flexibility and coordination. It is suitable for rehabilitation after hand injury.

## PACKAGING

- One Digi-Flex® of the selected colour / resistance
- Set of 5 different colours/resistances (with display)



Colour	Resistance: 1 digit / all digits	Reference code	Unit
	Yellow   0,7 kg / 2,3 kg	O2-O1O1O2	1 pc
	Red   1,4 kg / 4,5 kg	O2-O1O1O3	1 pc
	Green   2,3 kg / 7,3 kg	O2-O1O1O4	1 pc
	Blue   3,2 kg / 10,4 kg	O2-O1O1O5	1 pc
	Black   4,1 kg / 14,1 kg	O2-O1O1O6	1 pc
	Set of 5 resistances in 1 display	O2-O1O1I1	Set of 5 pcs
	Empty display	O2-O1O121	1 pc

# Digi-Extend®



The Digi-Extend® is a unique therapy device that can isolate the extension of all five fingers and each of the individual finger joints. The product strengthens the intrinsic and extrinsic muscle groups to maximise extensor tendon gliding, joint motion, coordination and muscle function in the hand and forearm.

Patients can choose between 5 resistance levels that are easily changeable. This way, you can vary the hand positions for specific exercises to help maximise your exercise goals and prevent injury.

Reference	Description	Unit
O2-O1O2O1	Digi-Extend   Bands (tan to green) included	1 pc
O2-O1O2I1	Digi-Extend Set   4 units   1 Display	set
O2-O1O2O2	Digi-Extend Band Set   25 pc 5x  5x  5x  5x  5x	set
O2-O1O221	Digi-Extend Empty Display (metal)	1 pc

Bands included:

- 5x
- 4x
- 3x
- 2x






## Handmaster Plus™

9 muscles close the hand, 9 muscles open the hand and 9 muscles control the position of the hand. The HandMaster Plus™ trains and balances all of these 27 muscles. Therefore, it is the ideal finger, thumb, hand and wrist exerciser!



SOFT | FOR LIGHT EXERCISE  
PURPLE BALL - YELLOW CORD



MEDIUM | FOR MEDIUM EXERCISE  
RED BALL - BLACK CORD



FIRM | FOR STRENGTH, SPEED TRAINING  
ORANGE BALL - GREY CORD

Reference	Description	Colour	Unit
O2-O5O1O2	Handmaster Plus   Purple ball - Yellow Cord		1 pc
O2-O5O1O3	Handmaster Plus Medium   Red ball - Black cord		1 pc
O2-O5O1O4	Handmaster Plus Firm   Orange ball - Grey cord		1 pc
O2-O5O1I1	Handmaster Plus   Set of 3 (one of each level)		Set

3-pack set  
=  
3 balls for the price of 2!

## Prohands® series

VIA™ is an entry level hand exerciser designed with young athletes in mind. It has a round, soft look and wide finger pistons for smoother function & improved finger performance.

Colour	Resistance (per digit)	Reference code	Unit
	Light   1,5 kg	O2-O6O1O2	1 pc
	Medium   2,5 kg	O2-O6O1O3	1 pc



### SPECIFICATIONS

Resistance: available in 2 different levels

PRO™ provides serious athletes with the most advanced hand exerciser in the world. The PRO™ features a larger profile, superior steel springs and an integrated ergonomic palm bar for comfort.

Colour	Resistance (per digit)	Reference code	Unit
	Medium   3,5 kg	O2-O6O3O3	1 pc
	Heavy   4,5 kg	O2-O6O3O4	1 pc
	Extra heavy   5,5 kg	O2-O6O3O5	1 pc



### SPECIFICATIONS

Resistance: available in 3 different levels



## VariGrip Therapy

This product has the ability to **adjust tension in each individual finger**. The resistance from each single finger knob can easily be adjusted incrementally on the patient's need, which minimises overcompensation and strain. This also allows the patient to progress steadily.

Comparison Chart (kg per digit)

VariGrip	X-Light	Light	Medium	Heavy	X-Heavy
Digi-Flex® / Flex-Ion	0,7	1,4	1,4 2,3	2,3 3,2	3,2 4,1

Each VariGrip's highest resistance ends where the next VariGrip resistance starts. Any other hand exerciser increases from light (1,4 kg) to medium (2,3 kg). That's almost a 70% increase gap, which can cause an overload on the muscles, resulting in injuries. The seamless progress from one VariGrip to the other, supported by the **padded finger and palm area**, guarantees a safe recovery and better feel for those with sensitivity issues.

Reference code	Colour	Resistance (range per digit) and level	Unit
O2-100102	Yellow	X-Light   0,7 - 1,4 kg   Level 1	1 pc
O2-100106	Black	X-Heavy   4,1 - 5,4 kg   Level 5	1 pc
O2-100107		VariGrip set   5 resistances   1 display	Set of 5
O2-100108		VariGrip empty display   Plastic	1 pc



## VariGrip Sport

The VariGrip Sport has the patented ability to **adjust the tension level in each individual finger**. This allows the user to create a custom programme to build and improve endurance, strength and dexterity more effectively for individual fingers, as well as the entire hand, wrist and forearm.

The VariGrip Sport is the only exerciser of its kind that includes a **callus builder**, ideal for climbers, gymnasts, gamers, musicians, and athletes of all levels.

This product also has **soft touch digits**, which provide comfort for users with sensitivity issues.

Reference code	Description	Resistance (kg)	Unit
O2-100303	VariGrip Sport	2,3 - 4,1 kg	1 pc



## Constant Force X-Tend

Unlike elastic-based exercisers that increase in resistance the further the unit is extended, this product uses retractable-based resistance, which keeps a constant level of resistance from start to finish. Ergonomic in design, the Constant Force X-Tend works the intrinsic and extrinsic extensors of the hand and can isolate MP/PIP/DIP joint motion.

Reference code	Colour	Resistance (per digit)	Unit
O2-100403	Red	Light   0,11 kg	1 pc
O2-100405	Blue	Heavy   0,27 kg	1 pc





## POWER-WEB


### Power-Web® Combo



The Power-Web® Combo features a unique advantage of two resistance levels in one, offering greater flexibility and value.

**The #1 hand exercise web in the world!**

**Each Power-Web® comes with a 39 exercise poster!**



#### SPECIFICATIONS

36 cm

450 gr

Reference code	Colour	Resistance	Unit
O2-O2O1O1		Tan - Red   Ultra light - Medium	1 pc
O2-O2O1O2		Yellow - Green   Light - Heavy	1 pc
O2-O2O1O3		Red - Blue   Medium - Super heavy	1 pc
O2-O2O1O4		Green - Black   Heavy - Ultra heavy	1 pc

## POWER-WEB

### Power-Web® Senior



The unique patented design allows a variety of exercises to be performed comfortably and effectively. Constructed of the finest rubber with special agents added for durability and strength, the Power-Web® will last for years.

#### SPECIFICATIONS

36 cm

450 gr

Reference code	Colour	Resistance	Unit
O2-O2O2O1		Tan   Ultra light	1 pc
O2-O2O2O2		Yellow   Light	1 pc
O2-O2O2O3		Red   Medium	1 pc
O2-O2O2O4		Green   Heavy	1 pc
O2-O2O2O5		Blue   Super heavy	1 pc
O2-O2O2O6		Black   Ultra heavy	1 pc



## POWER-WEB

### Power-Web® Junior



Just half the size of the Power-Web® Senior, this product has all the same features as its larger counterpart.

#### SPECIFICATIONS

18 cm

225 gr

Reference code	Colour	Resistance	Unit
O2-O2O3O2		Yellow   Light	1 pc
O2-O2O3O3		Red   Medium	1 pc
O2-O2O3O4		Green   Heavy	1 pc





## Power-Web® Flex-Grip MD CE

A distinctive and accessible exerciser that fits in the palm of your hand and provides the ability to work dexterity, extension and strengthening in your hands, fingers and thumbs. The Power-Web® Flex-Grip® has finger and thumb holes and has two sides: one side featuring a ball to squeeze and the other side featuring buttons and a thumb pad to exercise fingers separately.



### SPECIFICATIONS

 Resistance: available in 3 different levels

 Material: Latex-free

Colour	Resistance	Reference code	Unit
	Yellow   Light	O2-O211O2	1 pc
	Red   Medium	O2-O211O3	1 pc
	Green   Heavy	O2-O211O4	1 pc



## Adjustable hand grip

Train your hand strength with the Mambo Max Adjustable Hand Grip. This product comes in two different versions. A light (orange) version and a heavy (grey) version. Both can be adjusted to the tension you desire by turning the control knob.



LIGHT VERSION



HEAVY VERSION

Description	Resistance	Reference code	Unit
Mambo Max Adjustable Hand Grip   <b>Light</b>	5 to 20 kg   11 to 44 lb	O2-O7O1O2	1 pc
Mambo Max Adjustable Hand Grip   <b>Heavy</b>	10 to 40 kg   22 to 88 lb	O2-O7O1O4	1 pc



## Foam hand grip

Due to its foam handles, the Mambo Max Foam Hand Grip provides a firm grip for the entire workout session.



Description	Gripping power	Reference code	Unit
Mambo Max Foam Hand Grip	20 kg   44 lb	O2-O8O1O1	Pair



Squeeze balls & eggs




The Moves Squeeze Ball provides a progressive sequence of 5 colours/resistances. It strengthens hands, fingers, forearms and will work the muscles separately or in tandem. It's a small and accessible product that you can bring anywhere.

The Moves Squeeze Egg, also available in 5 levels, is ergonomically shaped to provide highly effective exercise for the rehabilitation and strengthening of weakened or injured fingers, hands and wrists.

Both the Balls and Eggs are compact and easy to use products that you can take everywhere you go. They can help to:

- ✓ Rehabilitate post-operative / post-traumatic hands and arms
- ✓ Restore function in hands, arms and shoulders following a stroke
- ✓ Reduce Pain and stiffness in arthritic fingers, hands and wrists
- ✓ Relax tense, sore muscles
- ✓ Reduce stress
- ✓ Strengthen ageing hands and arms




Colour	Resistance	Reference code - Ball	Reference code - Egg	Unit
	Yellow   Extra soft	02-030102	02-030202	1 pc
	Red   Soft	02-030103	02-030203	1 pc
	Green   Medium	02-030104	02-030204	1 pc
	Blue   Firm	02-030105	02-030205	1 pc
	Black   Extra firm	02-030106	02-030206	1 pc

PACKAGING

- ✓ 1 individually packed pc of the Moves Squeeze Ball or Egg

Wrist exerciser

The Moves Wrist Exerciser is an adjustable wrist and forearm trainer which allows the user to exercise with progressive resistance. Its ergonomically designed forearm support keeps the wrist in the right place. Together with the soft handgrip, adjustable to 4 positions, it ensures maximum comfort while exercising your wrist and forearm.



Description	Reference code	Unit
Moves Wrist Exerciser	02-070201	1 pc



## Comfort putty



The Moves Comfort Putty strengthens weakened grasp through exercise and is used for rehabilitation or stress relief.

It is odourless, does not stain and will not leave any residue on the patient's hands. Moreover, the Moves Comfort Putty is non-toxic, greaseless, hypo-allergenic and it contains **antimicrobial agents**.

Be assured of a top-tier product due to the Moves label!



### SPECIFICATIONS

- Resistance: available in 6 different levels
- Material: Latex-free, Gluten-free, contains anti-microbial agents
- Available in 2 weights: 85 g and 454 g

Colour	Resistance	Reference code		Unit
		85 g	454 g	
	White   Special Soft	O2-O9O1O1	O2-O9O4O1	1 pc
	Tan   Extra soft	O2-O9O1O2	O2-O9O4O2	1 pc
	Yellow   Soft	O2-O9O1O3	O2-O9O4O3	1 pc
	Red   Medium	O2-O9O1O4	O2-O9O4O4	1 pc
	Green   Firm	O2-O9O1O5	O2-O9O4O5	1 pc
	Blue   Extra firm	O2-O9O1O6	O2-O9O4O6	1 pc



### PACKAGING


- 1 putty of the selected colour/resistance pre-packed in a container of the chosen weight

## Why antimicrobial agents?

Antimicrobial agents **kill bacteria** without additives and they **prevent the growth of odour-causing microbes**. Moves Comfort Putty contains natural anti-microbial agents without the addition of toxic chemicals. This ensures that the Moves Comfort Putty material **retains its hygienic, non-greasy feel** from day one. Cheaper alternatives on the market do not contain these agents, which makes them the ideal breeding ground for bacteria and fungi. Why? Because a small number of bacteria and other micro-organisms will always be present on the patient's hands and because putty should be stored in room temperature conditions. This makes antimicrobial agents **a must for this product**. Apart from protecting the putty from bacteria and fungi, antimicrobial agents also **extend the product's lifespan** and **improve its cleanliness**.

### Benefits of antimicrobial agents:

- Non-greasy
- Hygienic
- Clean
- Durability
- Odourless
- Allows multiple users









New: Therapie Sand



MDCE

Thérapie Sand is an innovative material that provides a soothing tactile-kinesthetic experience and is equally popular among both children and adults. The sand promotes fine motor skills and sensitivity, stimulates circulation to relieve everyday stress, and restores mobility and overall strengthening of the wrists and ankles.

This product is primarily used in therapy to promote sensitivity and optimize the mobility of finger, hand, feet, toe, and ankle joints. It improves overall strength and hand-foot coordination. Additionally, the Thérapie Sand is an excellent educational aid, particularly in the fields of occupational therapy and physiotherapy, for the purpose of stimulating creativity and development.

With the Thérapie Sand you can train, play and destress all at once.

SPECIFICATIONS

-  Material: Latex-free sand
-  Available in 2 different weights: 454 g - 2 kg

Colour	Reference code	Description	Unit
	O2-24O4OI	Thérapie Sand   454 g   Grey	1 pc
	O2-24O5OI	Thérapie Sand   2 kg   Grey	1 pc





# O3 Upper & Lower Body







## Oxycycle 1 - active magnetic



For ACTIVE training only | For HOME use only

The Moves Oxycycle 1 Pedal Exerciser is used for upper and lower body exercises. This is a great product to progressively strengthen your arms and legs, because you can **adjust the resistance**.

### SPECIFICATIONS

- ↔ Measurements: 50 x 46 x 38 cm
- 📏 Weight: 9 kg
- 📄 Multi-function LCD display readout of scan, time, distance, calories and RPM
- ✔ Including: non-skid mat, interchangeable handgrip pedals and velcro strap pedals
- 🕒 **Maximum 1 hour per day**  
**Maximum 30 minutes per session**

Reference code	Description	Unit
O3-O1O1O1	Moves Oxycycle 1 Pedal Exerciser	1 pc



## Oxycycle 2 - passive



For PASSIVE training only | For HOME use only

The Moves Oxycycle 2 Pedal Exerciser is a great tool for personalised, low-stress workouts at home. The **adjustable speed level** makes it the ideal product to gradually strengthen your arms and legs. This pedal exerciser provides low-impact exercise to increase strength and flexibility. It's versatile for both upper and lower body and it improves the blood circulation.

It offers up to 60 Watts of motorized assistance, whereas cheaper models only carry a 30-Watt motor. So the Oxycycle 2 has enough power to maintain a constant rotational speed. A safety stop will avoid injury in case of cramp or muscle spasm.

### SPECIFICATIONS

- ↔ Measurements: 50 x 46 x 38 cm
- 📏 Weight: 8,5 kg
- 📄 Multi-function LCD display readout of scan, time, count, RPM-level and estimated calories burned  
**Hz: 50**  
**RPM: 25 - 65 (Rotations Per Minute)**
- ✔ Including: non-skid mat, interchangeable handgrip pedals and half-foot pedals
- 🔌 This product comes with a UK or EU plug
- 🕒 **Maximum 1 hour per day**  
**Maximum 30 minutes per session**

Reference code	Description	Unit
O3-O1O1O2	Moves Oxycycle 2 Pedal Exerciser   EU-plug	1 pc
O3-O1O1O2 UK	Moves Oxycycle 2 Pedal Exerciser   UK-plug	1 pc





## Oxycycle 3 - active & passive CE

For ACTIVE and PASSIVE training | For HOME use only

Increase your strength and flexibility with the Moves Oxycycle 3! This low-impact pedal exerciser can **adjust both the speed level and resistance**, which makes it the ultimate product to build strength in your arms and legs.

Unlike most cheaper models, which only carry a 30-Watt motor, the Moves Oxycycle 3 offers up to 80 Watts of motorized assistance. This assures a constant rotational speed. For your safety, a stop is integrated to prevent injury in case of a cramp or muscle spasm.

### SPECIFICATIONS

- ↔ Measurements: 50 x 46 x 38 cm
- 📏 Weight: 8,5 kg
- 📄 Multi-function LCD display readout of scan, time, count, RPM-level and estimated calories burned
  - Hz: 50
  - RPM: 25 - 65 (Rotations Per Minute)
- ✓ Including: non-skid mat, interchangeable handgrip pedals and full-foot pedals
- 🔌 This product comes with a UK or EU plug
- 🕒 **Maximum 1 hour per day**  
**Maximum 30 minutes per session**

Reference code	Description	Unit
O3-O1O1O3	Moves Oxycycle 3 Pedal Exerciser   EU-plug	1 pc
O3-O1O1O3 UK	Moves Oxycycle 3 Pedal Exerciser   UK-plug	1 pc



## Oxycycle pedals

All Moves Oxycycle Pedals are interchangeable and fit the Oxycycle 1, 2 and 3!

Reference code	Description	Unit
O3-O2O1O1	Moves Velcro Strap Pedals	pair
O3-O2O1O2	Moves Half-Foot Pedals	pair
O3-O2O1O3	Moves Full-Foot Pedals	pair
O3-O2O2O1	Moves Handgrip Pedals	pair



VELCRO STRAP PEDALS



HALF-FOOT PEDALS



FULL-FOOT PEDALS



HANDGRIP PEDALS




# CuraMotion

The CuraMotion is used for the improvement of muscle strength and joint range of motion.


The CuraMotion Exer Series has an adjustable resistance and a multi-functional LCD-display readout scan which indicates a lot of interesting metrics. All products are easily adjustable with a quick release knob.

## SPECIFICATIONS

-  Multi-functional LCD-display readout of scan, rotation count, rotations per minute, time exercised and calories burned.

Reference code	Description	Unit
O3-O3O1O1	CuraMotion Exer 1	1 pc
O3-O3O1O2	CuraMotion Exer 2	1 pc
O3-O3O1O3	CuraMotion Exer 3	1 pc
O3-O3O1O4	CuraMotion Exer 4	1 pc

### EXER 1


 Measurements:  
96 x 45 x 45 cm

 Weight: 11,5 kg

 Wall mount included



### EXER 2


 Measurements:  
146 x 45 x 54 cm

 Weight: 13,5 kg

 Wall mount included  
 Adjustable height



### EXER 3


 Measurements:  
37 x 38 x 70 cm

 Weight: 9 kg

 With foot stands



### EXER 4

 Measurements:  
96 x 45 x 45 cm

 Weight: 13 kg

 Wall mount included







### Foldable pedal exerciser MD CE

The Moves Pedal Exercisers are optimal for restoring impaired muscle strength, improving blood circulation and a better coordination for the legs and arms. It can be used on a table top as an arm exerciser or on the floor for a leg workout. A convenient resistance control knob allows a broad range of pedal resistance for a progressive exercise programme. This Moves Pedal Exerciser is foldable for compact shipping and storage. This product also provides an easy setup.

#### SPECIFICATIONS

- ↔ Measurements: 49 x 36 x 27 cm
- 📏 Weight: 2,3 kg
- 📄 Multi-function LCD display readout of scan, time, RPM, total rotations and calories
- ✔ Including: pedal straps that are adjustable to 2 positions

Reference code	Description	Unit
O3-OIO2O3	Moves Foldable Pedal Exerciser + LCD Display	1 pc



### Standard pedal exerciser | premium chrome MD CE

This Moves Pedal Exerciser is pre-assembled. The product has a high standard chrome finish and comes in a colour box.

#### SPECIFICATIONS

- ↔ Measurements: 52 x 42 x 26 cm
- 📏 Weight: 2 kg
- ✔ Option to buy 1 pc or 6 pcs in bulk
- ✔ Including: pedal straps that are adjustable to 2 positions

Reference code	Description	Unit
O3-OIO2O1	Moves Standard Pedal Exerciser   Premium Chrome	1 pc



### Standard pedal exerciser K/D MD CE

Are you looking for a budget solution to stay active after an injury? Or do you want a product to improve your coordination and the blood circulation of the upper and lower body? Then the Knock-down Moves Pedal Exerciser is perfect for you. It only takes 6 screws to (dis)assemble this product, so this product is ideal for compact shipping and storage.

#### SPECIFICATIONS

- ↔ Measurements: 51 x 41 x 28 cm
- 📏 Weight: 2 kg
- ✔ Including: pedal straps that are adjustable to 2 positions

Reference code	Description	Unit
O3-OIO2O2	Moves Standard Pedal Exerciser K/D   White	1 pc







## Fit stretch MD CE

The Fit Stretch is a fantastic, simple and portable tool to stretch your lower leg muscles. It is the ultimate product for warm-ups and cool-downs. The use of a Mambo Max Fit Stretch helps increase your ankle, shin and calf flexibility.

### SPECIFICATIONS

- 🔍 Material: made of durable moulded plastic
- ✅ Including non-skid foot pad



Reference code	Description	Unit
O3-O5OIOI	Mambo Max Fit Stretch   Single Leg	1 pc



## Wooden incline board MD CE

The Incline Board increases extension in the lower and upper back while stretching the calves. It can be used to:

- improve the general posture,
- stretch the calf muscles of athletes, football players and sportspersons in general,
- stretch the upper-leg muscles (Ischias),
- increase extension in the lower and upper back.

### THE MOVES WOODEN INCLINE BOARD - FIX

The perfect product if you are looking for an easy and budget-friendly Incline Board. It's made out of polymerised wood for a strong feel and high user weight. The inclination (17,5°) is perfect for physiotherapy exercises.

↔ Measurements: 45 x 30 x 16 cm

Reference code	Description	Unit
O3-O4O0O3	Moves Wooden Incline Board - Fix	1 pc







### New: Adjustable wooden incline board

The Mambo Max adjustable wooden incline board completes your workout setup. It is crafted from durable poplar wood and has a multi-layered board. This board can be adjusted in 4 different angles and conveniently folds for easy storage. Its slip-resistant stepping surface ensures stability during workouts. Experience enhanced stretching, rehabilitation, and strength training with this versatile incline board, perfect for targeting various muscle groups and improving flexibility and balance.

#### SPECIFICATIONS

- ↔ Measurements: 46 x 38 x 5,8 cm
- Q Material: poplar wood
- ∠ Inclination: adjustable to 10°, 20°, 30°, or 40° angles

Reference code	Description	Unit
O3-O4O0O2	Mambo Max adjustable wooden incline board	1 pc



### Plastic incline board MD CE

Apart from stretching the calf and upper-leg muscles, the incline board is a fantastic device to improve general posture and extension in the lower back. This Incline Board is made of a composite plastic for a maximum load of 250 kg. It adjusts to 4 levels (5°, 15°, 25° and 35°) and folds flat (8 cm) to store easily. The last slot on the board makes an excellent hand-hold for stretching. The non-slip top measures 35 x 35 cm.

Reference code	Description	Unit
O3-O4O1O2	Incline Board Standard	1 pc





# O4 Wellness





Mini massage gun €€

The Moves Mini Massage Gun is a small sized, lightweight, silent and cordless massager that can be used anywhere, anytime. Depending on your preference of frequency and force, you can adjust the intensity of your massage using the 4 speed levels. The Moves Mini Massage Gun provides relief for muscle stiffness and soreness. Furthermore it reduces tenderness in neck, shoulders, waist and legs, while also accelerating your blood circulation. This product comes with 5 attachments that can be switched to target specific problem areas, such as your neck, shoulders, back, legs, arms, hands and feet.

SPECIFICATIONS

↔ Measurements: 14,6 x 9,8 x 4,6 cm

🛒 Weight: 626 g



- **Speed levels:** 4
- **Motor type:** brushless motor
- **Motor:** 31,2 W / 50 dB
- **Rated speed:** 1800 - 3200 BPM
- **Amplitude:** 10 mm
- **Pressure strength:** 15 kg



- **USB-C charging cable**
- **Battery capacity:** 2600 MAh.
- **Charger power output:** 5V2A
- **Charge time:** 3 - 4 hours
- **Use after fully charging:** up to 4 hours



Reference code	Description	Unit
04-080101	Moves Mini Massage Gun	1 pc



1. SPHERICAL HEAD



2. BULLET HEAD



3. U-SHAPED HEAD



4. SPADE HEAD



5. PLUM BLOSSOM HEAD

1. **Spherical head** - suitable for large muscle groups such as arms, waist, hips, thighs and calves
2. **Bullet head** - suitable for deep tissues such as acupoints, tendons, joints, palms and soles of feet
3. **U-shaped head** - suitable for both sides of the spine and Achilles tendon
4. **Spade head** - suitable for relaxation and shaping of various muscle parts
5. **Plum blossom head** - suitable for deeper parts of the legs, buttocks and large areas of muscle

PACKAGING

- ✓ Cardboard packaging
- ✓ Including a convenient storage case, USB-C charging cable and 5 attachments





## Triggerpoint rollers

Massage rollers are used to put pressure on muscles' triggerpoints and to loosen the tension in the muscles. Triggerpoints are often the cause of different pains.

### 1. SPIKY FOAM ROLLER


This product is designed for **deeper massage & acupressure** therapy. The spikes make it possible to get deeper into the muscles to relieve tension and soreness and to improve mobility. The Spiky Foam Roller has a bi-directional zone for comfort when rolling along the spine. It is mainly used for stretching, self-massage, body support and stability exercises.

#### SPECIFICATIONS

 Length: 33 cm

 Diameter: 13 cm

Reference code	Description	Unit
O4-O5O1O3	Mambo Max Spiky Foam Roller	1 pc



### 2. HOLLOW FOAM ROLLER

The Hollow Foam Roller easily helps people who are experiencing **long term pains**. Its 3-dimensional pattern with different shapes replicates the massage from a therapist's hands. The larger surfaces allow a superficial massage while the smaller shapes are meant for deeper massage.

#### SPECIFICATIONS

 Length: 33 cm

 Diameter: 14 cm

Reference code	Description	Unit
O4-O5O1O1	Mambo Max Hollow Foam Roller	1 pc



### 3. ERGONOMIC FOAM ROLLER

The Mambo Max **Ergonomic Foam Roller** is a curved, hollow foam roller, which enables a lot of exercises. Besides rolling back and forth, it is possible to make sideways movements. More muscles will be addressed through which the body will be toned even more. The soft foam in combination with the nipples is perfect for trigger point massages. The central groove will **reduce pressure on the spine**.

#### SPECIFICATIONS

 Length: 33 cm

 Diameter: 15 cm

Reference code	Description	Unit
O4-O5O1O2	Mambo Max Ergonomic Foam Roller	1 pc







4. 3-IN-1 FOAM ROLLER

The Mambo Max 3-in-1 Foam Roller combines 3 different foam rollers in one package. Thanks to its smooth surface, the **outer part** is perfect for a superficial massage. The **second part** has grooves which can be used for a deeper massage and to exercise higher pressure on the trigger points. Finally, with the **inner massage stick** it's possible to improve blood circulation and the muscles' recovery time.

This last item, the **Mambo Max Massage Stick**, is also available as a separate product

SPECIFICATIONS

-  Length: all 3 parts are 46 cm
-  Diameter:
- outer part: 15 cm
  - second part: 10 cm
  - inner massage stick: 4 cm

Reference code	Description	Unit
O4-O5O3O1	Mambo Max 3-in-1 Foam Roller	1 set
O4-O5O2O1	Mambo Max Massage Stick	1 pc



MASSAGE STICK





**IASTM tools**

Moves IASTM Tools (Instrument-Assisted Soft Tissue Mobilization Tools) enable therapists to **treat acute and chronic soft tissue dysfunction**. The tools are used in physical therapy, trigger point therapy, for myofascial release and various other treatments. They can be used on many body parts, including muscles, tendons, ligaments, and fascia.



IASTM works like a standard massage, but specialized instruments are used instead of a physical therapist's hands. These instruments are designed to apply a specific amount of pressure to tissue so that a therapist can work through it without causing damage or pain to the client.

**HOW TO CHOOSE THE RIGHT IASTM TOOL?**

Each tool has multiple shapes and finished edges to increase its employability. The smaller ones are generally used for the treatment of smaller parts of the body, smaller strokes and reaching superficial tissue structures. The larger and heavier ones are used to reach deeper and are suitable for the treatment of larger parts of the body. Some instruments are also used to “feel” blockages by stroking/scanning across the tissue. The tips of these instruments can be either blunt or sharp, depending on the treatment given. As a general rule, the sharper the tool, the more intense and specific the soft tissue treatment. Finally, the choice of the right tool is different for every therapist. Everyone has their personal preferences and the tool should feel comfortable in one's hand. With our 6 tools, you can easily customize your treatment based on both your needs and the clients needs.

**SPECIFICATIONS**

- ✓ Sold as a set of 6 (different types)
- Q Material: high-grade stainless steel and a brushed finish (less prone to scratches and smudges are less visible)



**USE AND DIMENSIONS OF EACH TOOL**

The first tool is ideal to work some of the **larger muscles**.  
E.g. the lower leg, chest, upper arms, back, lats and upper- or lower legs.  
20 x 3 cm

The second instrument is often used for **medium sized areas**.  
E.g. the neck, calves or pectoral muscles.  
19 x 4 cm

The largest tool has two handles and a round side for deep treatment. It is designed for **large muscle groups**.  
E.g. the back, upper arms, hamstrings and thighs.  
45 x 4 cm

The fourth instrument is small, yet versatile. It's mostly used for working with **small muscles**.  
E.g. the treatment of arms, shoulders, hands, feet and ankles.  
14 x 3 cm

This sturdy and versatile tool is perfect for **reducing tension, muscle pain and rebalancing the fascial network and body**. This tool is also useful for **scanning tissue**.  
7 x 6 cm

Finally, with this easy to hold tool you can treat **various muscle groups and joints**.  
E.g. the hands, forearms, chest, back, upper- and lower leg.  
9 x 7 cm

Reference code	Description	Unit
O4-O9O111	Moves IASTM Tools	Set of 6



## Massage tools

### THE KNOBBLE® II

This is a massage tool made of durable polymer with a non-latex band. The Knobble® II helps you to put more pressure on the area that you want to massage. Use either the tip or the base for applying pressure on yourself or others.



### THE JACKNOBBER® II

It's a massage tool that can be used in many ways to apply deep pressure to myofascial triggerpoints and areas of muscle tension. The smaller knobs provide more intense pressure, while the larger knobs give you a more diffuse pressure. Use one knob at a time to press deeply into taut bands with embedded triggerpoints. Use two knobs at a time on either side of the spine.




### THE INDEX KNOBBER® II

The Index Knobber® II is a simple device designed to let you apply deep pressure to muscles within easy reach of your hand. It is the ideal instrument for clinic or home use allowing for firm, sustained pressure on targeted muscles. Its unique design allows it to be used in several hand positions and it can easily be used with both your right and left hand.



### THE PALMASSAGER™

This tool is one of the most efficient massage therapy aids available anywhere. Hold it firmly in the palm of your hand and use one, two or three of its round knobs to massage sore muscles. It is used for back, shoulders, arms and legs.




### THE BACKNOBBER® II

The Backnobber® II massage tool allows the user to apply as much pain-releasing pressure as they require. Remarkable effects can be achieved using this device. A locking mechanism permits it to come apart for convenient traveling and storage.



Reference code	Description	Unit
O4-O311O1	Jacknobber II	1 pc
O4-O311O2	Index Knobber II	1 pc
O4-O311O3	PalmaSSager	1 pc
O4-O311O4	Backnobber II	1 pc
O4-O311O5	Knobble II	1 pc







## Massage ball, roll & ring

These massage tools are used for relaxing muscles and reducing tension. The Ball, Roll and Ring help to strengthen hands and feet and to improve coordination and flexibility. They can be used for massage as well as to improve dexterity skills.

Reference code	Description	Unit
O4-O3O1O3	Mambo Max Massage Ball   6 cm	1 pc
O4-O3O1O5	Mambo Max Massage Ball   7 cm	1 pc
O4-O3O1O7	Mambo Max Massage Ball   8 cm	1 pc
O4-O3O1O9	Mambo Max Massage Ball   9 cm	1 pc
O4-O3O1I1	Mambo Max Massage Ball   10 cm	1 pc
O4-O3O1I2	Mambo Max Massage Roll   16 cm	1 pc
O4-O3O1I3	Mambo Max Massage Ring	1 pc



## moves

### Sport tape



The Moves Sport Tape is a non-elastic tape used **to limit motion**. The tight binding nature of the Sport Tape **decreases circulation**. This tape is made from quality cotton fibres which makes it soft and comfortable, while still providing maximum compression. The skin-friendly Moves Sport Tape is latex-free, breathable and water-resistant. The special adhesive (hypoallergenic Zinc oxide hot melt adhesive) makes sure that the tape is fixed and secure during exercise. The serrated edges make it very easy for the tape to be torn by hand.

#### SPECIFICATIONS

Dimensions: 3,8 cm x 10 m

Available in white

Reference code	Description	Unit
O4-O7O1O1	Moves Sport Tape	1 pc



## moves

### Kinesiology tape



Moves Kinesiology Tape is used **to facilitate motion, inhibit pain and improve circulation** while still supporting the muscles and joints. With its wave-like structure, this tape can be stretched up to 170% lengthways. This tape is also made from fine cotton fibres, which results in a soft and comfortable tape. The skin-friendly product is hypoallergenic, latex-free, breathable and water-resistant. The adhesive side is made of acrylic which prevents sweat and moisture from penetrating the tape.

#### SPECIFICATIONS

Dimensions:

- 5 cm x 5m (in all colours)
- **New:** 5 cm x 35m (in tan)

Available in 4 colours (tan | black | pink | light blue)

Colour	Description	Dimensions	Reference code	Unit
	Moves Kinesiology Tape   <b>Tan</b>	5 cm x 5 m	O4-O6O1O1	1 pc
	Moves Kinesiology Tape   <b>Black</b>	5 cm x 5 m	O4-O6O1O2	1 pc
	Moves Kinesiology Tape   <b>Pink</b>	5 cm x 5 m	O4-O6O1O3	1 pc
	Moves Kinesiology Tape   <b>Light blue</b>	5 cm x 5 m	O4-O6O1O4	1 pc
	Moves Kinesiology Tape   <b>Tan</b>	5 cm x 35 m	O4-O6O1O5	1 pc





## Sitting cushions

Mambo Max Sitting Cushions are air-filled cushions for active and dynamic sitting, recommended by therapists. They're easy-to-carry cushions and they activate intrinsic trunk muscles to support the spine. They relieve fatigue and adjust your posture. The surface from these sitting cushions can easily be cleaned with soap and water. All cushions can be used **as a seat** as well as **a balance trainer**.

### COMFORT CUSHION 33

Soft and comfortable surface

- ↕ Height: 7 cm
- ∅ Diameter: 14 cm
- ✓ Packed in a colour box, including Needle Pump
- ✓ Packed in a polybag only, doesn't include Needle Pump




NEEDLE PUMP

Reference code	Description	Unit	Icon
O4-O2O3O1	Mambo Max Comfort Cushion 33   Blue   + Needle Pump	1 pc	
O4-O2O3O2	Mambo Max Comfort Cushion 33   Red   + Needle Pump	1 pc	
O4-O2O3O3	Mambo Max Comfort Cushion 33   Black   + Needle Pump	1 pc	
O4-O2I3O1	Mambo Max Comfort Cushion 33   Blue   In polybag	1 pc	
O4-O2I3O2	Mambo Max Comfort Cushion 33   Red   In polybag	1 pc	
O4-O2I3O3	Mambo Max Comfort Cushion 33   Black   In polybag	1 pc	

### STANDARD CUSHION 37

One smooth side | One senso side (with knobs)

- ↕ Height: 7 cm
- ∅ Diameter: 38cm
- ✓ Including Faster Blaster pump



FASTER BLASTER PUMP

Reference code	Description	Unit	Icon
O4-O2O1O3	Mambo Max Standard Cushion 37   + Faster Blaster Pump	1 pc	

### STANDARD CUSHION 33

One smooth side | One senso side (with knobs)

- ↕ Height: 7 cm
- ∅ Diameter: 33cm
- ✓ Including Needle Pump



SMOOTH SIDE

Reference code	Description	Unit	Icon
O4-O2O1O1	Mambo Max Standard Cushion 33   + Needle Pump	1 pc	



**SENSAIR CUSHION**

One smooth side | One senso side (with knobs)

- ↕ Height: 7 cm
- ∅ Diameter: 31cm
- ✓ Including Needle Pump



SMOOTH SIDE


Reference code	Description	Unit
O4-O2O1O2	Mambo Max SensAir Cushion   + Needle Pump	1 pc



**SITTING WEDGE**

One smooth side | One senso side (with knobs)

- ↕ Height: 7 cm
- Width: 34cm
- ✓ Including Faster Blaster Pump



SMOOTH SIDE

Reference code	Description	Unit
O4-O2O2O1	Mambo Max Sitting Wedge   + Faster Blaster Pump	1 pc






## New: Adjustable ball chair

The Mambo Max adjustable ball chair offers you the ideal way to upgrade your work setup. This chair enables active sitting throughout the day. That way, you put an end to sedentary behavior in the workplace or at home. The ergonomic back support assists with spine alignment, it strengthens the core, and it corrects your sitting posture.

### SPECIFICATIONS

- Instability-induced motion helps to firm your core, improve circulation and posture
- Chair height ranges from 51 cm to 58 cm (20" to 23") to fit most standard desks
- Stable, 360° swivel base allows for easy movement
- Lightweight and compact
- Lockable easy glide caster wheels
- Anti-burst, removable oval-shaped ball seat, different from other ball chairs
- Max weight: 127 kg (280 lbs)

Colour	Description	Reference code	Unit
	Mambo Max Adjustable Ball Chair Black / Grey ball	O4-O2O4O1	1 pc



## Pilates

### FOAM ROLLER

The Mambo Max Foam Roller is used for Pilates training and spine exercises, more specifically back aches. It's ideal for balance and stabilisation exercises. Foam rollers are used for stretch and strength exercises, as well as for massage and relaxation.

### SPECIFICATIONS

- ↔ Length: 45 cm | 90 cm
- ∅ Diameter: 15 cm
- 🔍 Material: EVA foam (soft, yet durable polymer of ethylene and vinyl acetate)

Reference code	Description	Size	Unit
O4-O1O1O2	Mambo Max Pilates Foam Roller   EVA	90 cm	1 pc
O4-O1O1O4	Mambo Max Pilates Foam Roller   EVA	45 cm	1 pc





## SOFT-OVER-BALL

The Mambo Max Soft-Over-Ball is an inflatable ball used for resistance work in Pilates. It is a wonderful tool for adding focus and energy to your workouts. In addition, when slightly deflated, it can be used as a soft prop under the head, neck, feet, back...  
Now available in 3 different colours and sizes, giving you a range of 9 Soft-Over-Balls

Reference code	Description	Size	Unit
O4-O1O1O5	Mambo Max Pilates Soft-Over-Ball   <b>Green</b>	17 - 19 cm	1 pc
O4-O1O115		21 - 23 cm	1 pc
O4-O1O125		25 - 27 cm	1 pc



Reference code	Description	Size	Unit
O4-O1O113	Mambo Max Pilates Soft-Over-Ball   <b>Blue</b>	17 - 19 cm	1 pc
O4-O1O123		21 - 23 cm	1 pc
O4-O1O1O3		25 - 27 cm	1 pc



Reference code	Description	Size	Unit
O4-O1O1O7	Mambo Max Pilates Soft-Over-Ball   <b>Black</b>	17 - 19 cm	1 pc
O4-O1O117		21 - 23 cm	1 pc
O4-O1O127		25 - 27 cm	1 pc



## RING

The Mambo Max Pilates Ring is a resistive ring with foam padded grips. This product is used for exercises focusing on strength, balance, flexibility and posture. The ring adds light resistance to Pilates exercises and helps to target specific muscle groups. It's often used to help tone the arms and thighs, and to increase the core strength challenge of an exercise

## SPECIFICATIONS

- Ø Diameter: 38 cm
- Q Material: foam covered fiberglass

Reference code	Description	Unit
O4-O1O1O1	Mambo Max Pilates ring	1 pc





## Stepping stones

The non-slip Mambo Max Stepping Stone is an amazing tool to freshen up and enhance your workout. It can be used by both children and adults for indoor and outdoor activities to improve dynamic balance, body awareness and coordination. These Stepping Stones can also serve as a foot massage tool while you're seated.

### SPECIFICATIONS

- Ø Diameter: 16 cm
- ✓ Vivid colours: yellow, orange, red, purple, blue and green
- ✓ 6 pcs of a different colour packed in a cardboard box

Reference code	Description	Unit
O4-O3O141	Mambo Max Stepping Stone	6-pack





## Yoga

Mambo Max offers different Yoga items. These can be bought separately or in a Yoga Set. The set consists of the following products:

- 1 PVC yoga mat | 173 x 61 x 0,4 cm
- 2 EVA foam yoga blocks | 7,5 x 15 x 23 cm
- 1 yoga strap with buckle | 180 x 3,8 cm

Reference code	Description	Unit
O4-O1O201	Mambo Max Yoga Mat	1 pc
O4-O1O202	Mambo Max Yoga Block	5-pack
O4-O1O203	Mambo Max Yoga Strap	1 pc
O4-O1O211	Mambo Max Yoga Set   1 Mat   2 Blocks   1 strap	Set





# O5 Fitness





## Gym mat

Mambo Max Gym Mats provide a soft touch and comfortable exercise environment. They are perfect for gym & fitness. They have a non-slip, well-cushioned surface with a unique closed cell structure. The mats are equipped with plastic rings that are 40 cm apart so they can be suspended while not in use. Mambo Max also offers a Gym Mat Wall Hanger that can hold up to 10 mats.

The Mambo Max Gym Mats are available in 3 different sizes and thicknesses. They are **latex-free** and can be used in water. Each mat comes with a handy **storage/carry strap**.

### COMFORT GYM MAT

measurements: 140 x 60 x 1 cm


colours: and



### XTRA COMFORT GYM MAT

measurements: 180 x 60 x 1,5 cm


colours: and and and




### XTRA COMFORT GYM MAT

measurements: 180 x 100 x 1,5 cm

colours: and and



Reference code	Description	Unit
O5-O30001	Mambo Max Gym Mat Wall Hanger   For 10 mats	1 pc
O5-O30303	Mambo Max Comfort Gym Mat   140 x 60 x 1 cm      Rings   Strap	1 pc
O5-O30304	Mambo Max Comfort Gym Mat   140 x 60 x 1 cm      Rings   Strap	1 pc
O5-O30101	Mambo Max Xtra Comfort Gym Mat   180 x 60 x 1,5 cm      Rings   Strap	1 pc
O5-O30102	Mambo Max Xtra Comfort Gym Mat   180 x 60 x 1,5 cm      Rings   Strap	1 pc
O5-O30103	Mambo Max Xtra Comfort Gym Mat   180 x 60 x 1,5 cm      Rings   Strap	1 pc
O5-O30104	Mambo Max Xtra Comfort Gym Mat   180 x 60 x 1,5 cm      Rings   Strap	1 pc
O5-O30201	Mambo Max Xtra Comfort Gym Mat   180 x 100 x 1,5 cm      Rings   Strap	1 pc
O5-O30202	Mambo Max Xtra Comfort Gym Mat   180 x 100 x 1,5 cm      Rings   Strap	1 pc
O5-O30203	Mambo Max Xtra Comfort Gym Mat   180 x 100 x 1,5 cm      Rings   Strap	1 pc





## Aerobic step

Step aerobics is a form of aerobics that is practiced worldwide by millions of people. The reason for its popularity is simple, it offers a great cardio workout without the need for extensive space or equipment.

Furthermore, it's low impact and it doesn't stress the joints as much as running or jogging.

The height of the Mambo Max Aerobic Step is adaptable, which allows you to **adjust the intensity** of the work-out according to your preferences.

### SPECIFICATIONS



Measurements: 96 x 33 cm



Height can be altered to **15 / 20 / 25 cm** by using blocks underneath the step



100 kg

Reference code	Description	Unit
05-070101	Mambo Max Aerobic Step	1 pc




15 CM



20 CM



25 CM



## Fabric elastic band

The Mambo Max Fabric Elastic Band consists of 10 consecutive loops. This great toning and workout tool is suitable for every kind of athlete, whether you're a beginner or a professional. There are 3 different resistances available, while the 10 consecutive loops allow you to spice things up as well. Intensify your physical training by using a lower amount of consecutive loops resulting in a higher resistance.

This product is durable and extremely versatile. You can use it for stretching, toning or even as an alternative for dumbbells as a weight training accessory.

### SPECIFICATIONS



Total length: 120 cm  
1 loop length: 11,5 cm



Resistance: available in 3 different levels

Colour	Resistance	Reference code	Unit
Red	Red   Light   7 kg	05-120107	1 pc
Green	Green   Medium   10 kg	05-120110	1 pc
Blue	Blue   Heavy   15 kg	05-120115	1 pc





## Balance board

The Mambo Max Balance Board is an effective tool for balance training and fitness exercises. Stand on the rounded platform and start balancing. It is also a great tool to strengthen your ankles and thigh muscles. For a more severe balance training, add the included screw-on top to your Mambo Max Balance Board.

### SPECIFICATIONS

Ø Diameter: 42 cm

↕ Height: 5 - 8 cm



Reference code	Description	Unit
O5-O4O2O1	Mambo Max Balance Board	1 pc



## Wooden balance boards

These Balance Boards are the ideal solution if you want to improve your balance on an **ecologic, high-end product**. They are often used for proprioception training, which challenges and develops proprioceptive receptors. Proprioception helps to stabilise a joint during static and dynamic functional tasks. Proprioception exercises can decrease your risk of injury.

The Wooden Balance Board is available in 3 versions: Uni-Directional, Standard & Sport. Thanks to their different form and balance point, these products have a very different character and level of difficulty.

### UNI-DIRECTIONAL

This board offers a **more predictable** and **confident training** as it can only pivot in 1 direction. Making it ideal for children and elderly to improve body stability without the fear of falling. When turning this board upside-down, the arches can be used as push-up grips.

□ Measurements: 40 x 40 cm

↕ Height: 10 cm



### STANDARD

The Balance Board Standard is lower than the Sport version. In addition, with its **smoother pivot point**, it reacts completely different than the Sport. It's a great middle step between the Uni-Directional and Sport Balance Boards.

□ Measurements: 40 x 40 cm

↕ Height: 6 cm

△ Max inclination: 12,5°



Reference code	Description	Unit
O5-O4OOO1	Mambo Max Wooden Balance Board - Standard	1 pc
O5-O4OOO3	Mambo Max Wooden Balance Board - Uni-directional	1 pc



## Balance board freeman

The Mambo Max Balance Board Freeman is ideal for balance training. Stand on the rounded platform and start balancing. This is a great way to strengthen your ankles and thigh muscles. The handgrips can be used for upper body balancing and workouts.

### SPECIFICATIONS

Ø Diameter: 40 cm

↕ Height: 8 cm



Reference code	Description	Unit
O5-O4O2O2	Mambo Max Balance Board Freeman	1 pc



## Body twister

If you're looking for a tool to **slim your thighs, waist and hips**, the Mambo Max Body Twister is the answer! While twisting your body, the disc below your feet turns. This trains your glutes, abdominals and legs. Continue to twist your body from side to side to tone all the muscles in your body. Furthermore, the Mambo Max Body Twister is an effective product to **enhance balance and coordination**.

### SPECIFICATIONS

Ø Diameter: 25 cm

↕ Height: 3 cm



Reference code	Description	Unit
O5-100101	Mambo Max Body Twister	1 pc




## Ekilibre

The Mambo Max Balance Plate "Ekilibre" is designed for balance training. **Combining balance pad with balance board**, it offers the best of both worlds. It improves your sense of equilibrium, coordination, reaction and motor skills. It can be used for lying, kneeling, sitting and stepping movement and it can be incorporated with other exercise equipment.

Reference code	Description	Unit
O5-O4O5O1	Mambo Max Ekilibre   Beginner	1 pc
O5-O4O5O2	Mambo Max Ekilibre   Intermediate	1 pc


### BEGINNER



□ Measurements:  
60 x 39 cm  
↕ Height: 9 cm



### INTERMEDIATE



Ø Diameter: 39 cm  
↕ Height: 9 cm





## Balance trainer

The Mambo Max Balance Trainer is used to improve your sense of equilibrium, coordination, reaction and motor skills. It is also ideal for muscle exercise. You can adjust the Balance Trainer by inflating and deflating the air through the valve. You can use it for kneeling, sitting and stepping movement. To spice up your rehabilitation, incorporate this product with other exercise equipment.



- Ø Diameter: 45 cm
- ↕ Height: 16 cm
- 🔍 2 equal sides
- 🧑 100 kg



- Ø Diameter: 60 cm
- ↕ Height: 14 cm
- 🔍 1 smooth side  
1 senso side
- 🧑 100 kg

Reference code	Description	Unit
O5-O4O1O2	Mambo Max Balance Trainer   45 cm	1 pc
O5-O4O1O3	Mambo Max Balance Trainer   60 cm	1 pc

## Balance pad

This **EVA-foam** Balance pad provides the comfortable exercise environment needed to perform reinforcing balance exercises. Its soft and consequently unstable surface challenges the user to work the foot-, ankle-, and leg muscles.

### AVAILABLE IN 3 VERSIONS:



- Rectangular shaped pad.
- 📏 Measurements: 50 x 40 x 6 cm
- Colour: Blue



- Rectangular shaped pad.
- 📏 Measurements: 50 x 40 x 6 cm
- Colour: Anthracite



- 2 Oval shaped pads, sold per pair.
- 📏 Measurements: 36 x 22 x 6 cm
- Colour: Blue

Reference code	Description	Unit
O5-O4O3O1	Mambo Max Balance Pad   Blue	1 pc
O5-O4O3O3	Mambo Max Balance Pad   Anthracite	1 pc
O5-O4O3O2	Mambo Max Balance Pad   Oval   Blue	Pair



## Balance pad pro

This Balance Pad is a great training tool to improve your balance, stability, flexibility, strength and motor skills. The Balance Pad Pro **creates more instability** than our regular one, so you constantly need to correct body posture, making the exercises just that bit tougher. This Balance Pad is made out of **PU-material** instead of EVA-foam (which most balance pads are made from), which makes it **more durable** and less prone to small damages and wear.

### SPECIFICATIONS

Measurements: 40 x 50 x 6 cm

Reference code	Description	Unit
O5-O4O3O4	Mambo Max Balance Pad Pro   Rectangular   Blue	1 pc



## AB gym ball

The Mambo Max **AB (Anti-Burst)** gym ball increases flexibility, improves strength, fitness and posture. Exercising with this ball also enhances your balance and coordination. All AB gym balls are Anti-Burst: if the ball gets punctured by a sharp object, the ball won't burst, thus reducing the risk of an injury.

### SPECIFICATIONS

Available in 6 different sizes

Material: ecofriendly PVC  
(Heavy metal content in compliance with EN 71 Standards.  
Phthalates content meets EU-directive 2005/84/EC)

200 kg

### PACKAGING

- 1 AB gym ball in a cardboard retail pack
- Includes:  
Faster Blaster Pump, Exercise Chart, Plug and Plug Remover
- Value pack: a polybag containing only the ball and plug



Colour	Ball Size	Body Height	Mambo Max AB gym ball Box + Pump	Mambo Max AB gym ball Polybag	Unit
	45 cm	< 125 cm	O5-O1O1O2	n/a	1 pc
	55 cm	125 - 150 cm	O5-O1O1O3	O5-O1O2O3	1 pc
	65 cm	150 - 175 cm	O5-O1O1O4	O5-O1O2O4	1 pc
	75 cm	175 - 185 cm	O5-O1O1O5	O5-O1O2O5	1 pc
	85 cm	185 - 200 cm	O5-O1O1O6	n/a	1 pc
	95 cm	> 200 cm	O5-O1O1O7	n/a	1 pc



## AB peanut ball

The Mambo Max **AB (Anti-Burst)** Peanut Ball helps **users who struggle with the instability** of a regular gym ball. The peanut shape offers extra sideways stability, which gives more confidence to start working out with gym balls.

The Mambo Max AB Peanut Ball offers a more predictable therapy or fitness training as it can only roll in one direction. This certainly is a huge benefit for children, elderly or people who suffered an injury. In addition, peanut balls are very popular amongst pregnant women and are very helpful devices during labour.

### PACKAGING

- ✓ 1 AB Peanut Ball in a cardboard retail pack
- ✓ Includes: Faster Blaster Pump and plug.



150 kg

Colour	Peanut Ball Size	Mambo Max AB Peanut Ball + Pump	Unit
<span style="color: red;">●</span>	50 x 100 cm	O5-O11103	1 pc



## Dynadome

Half sphere, inflatable exercise dome. This versatile product enables the user to perform exercises on both sides of the product, allowing it to be used for a full body workout as well as for balance exercises. Basically, it's Mambo Max's economical answer to the expensive balance trainers on the market, but **with extra tubing and handles!**


### SPECIFICATIONS

- Ø Diameter: 59 cm
- ↕ Height: 21 - 25 cm
- ✓ Includes: **Faster Blaster Pump**
- ✓ Includes: **Tubing with handles**



175 kg

Reference code	Description	Unit
O5-O40401	Mambo Max Dynadome   Tubing with handles   Faster Blaster Pump	1 pc





## Ball accessories



ELECTRIC COMPRESSOR PUMP



HAND PUMP




FASTER BLASTER PUMP



NEEDLE PUMP



FOOT PUMP



BALL STACKER  
WITH VELCRO  
CLOSURE




BALL STABILISER THAT FACILITATES BALANCE.  
ITS NOTCHES ENABLE COMBINED AB GYM BALL / MOVES  
BAND EXERCISES.

Reference code	Description	Unit
O5-O2O1O1	Mambo Max Ball Stabiliser	1 pc
O5-O2O1O2	Mambo Max Ball Stacker	3 pcs
O5-O2O2O1	Mambo Max Faster Blaster Pump   2 x 120 cc	1 pc
O5-O2O2O2	Mambo Max Hand Pump   2 x 2000 cc	1 pc
O5-O2O2O3	Mambo Max Needle Pump   100 cc	1 pc
O5-O2O2O4	Mambo Max Foot pump	1 pc
O5-O2O2O5	Electric Compressor Pump	1 pc



# O6 Strength






## SoftMeds

Mambo Max SoftMeds are colour coded hand weighted balls that are small and flexible enough to grasp with one hand, yet big enough to hold with both hands. The soft, pliable vinyl cover makes it easy to hold and is a safe and effective way to add weight training to your workout.

Mambo Max SoftMeds are fully inflated and ready to use, straight out of the box! When needed, air pressure can be adjusted with a needle pump.



### SPECIFICATIONS

Ø Diameter: 12 cm

Reference code	Weight	Description	Colour	Unit
O6-O4OIO2	0,5 kg	Mambo Max SoftMeds   Tan		1 pc
O6-O4OIO3	1,0 kg	Mambo Max SoftMeds   Yellow		1 pc
O6-O4OIO4	1,5 kg	Mambo Max SoftMeds   Red		1 pc
O6-O4OIO5	2,0 kg	Mambo Max SoftMeds   Green		1 pc

## Coated dumbbell


Mambo Max Dumbbells are cast iron weights, coated with vinyl or neoprene. The vinyl coating is easy to wash or can be treated with a water based disinfectant. The neoprene coating has an added resin to give them a scrubby feel and superior grip.

### SPECIFICATIONS

- ✓ Available in 7 different weights + matching colour code
- ✓ The Dumbbell Rack can hold 6 Dumbbells



NEOPRENE DUMBBELLS



RACK WITH VINYL DUMBBELLS

Description	Colour	Weight	Reference code	Reference code	Unit
			Vinyl	Neoprene	
Mambo Max Dumbbell   Pink		0,5 kg	O6-OIOIO2	O6-OIO2O2	pair
Mambo Max Dumbbell   Yellow		1,0 kg	O6-OIOIO3	O6-OIO2O3	pair
Mambo Max Dumbbell   Red		1,5 kg	O6-OIOIO4	O6-OIO2O4	pair
Mambo Max Dumbbell   Green		2,0 kg	O6-OIOIO5	O6-OIO2O5	pair
Mambo Max Dumbbell   Blue		3,0 kg	O6-OIOIO7	O6-OIO2O7	pair
Mambo Max Dumbbell   Black		4,0 kg	O6-OIOIO8	O6-OIO2O8	pair
Mambo Max Dumbbell   Grey		5,0 kg	O6-OIOIO9	O6-OIO2O9	pair
Mambo Max Dumbbell Rack (empty)			O6-000IO1		1 pc



## Medicine ball

A Medicine Ball (also known as an exercise ball, a med ball, or a fitness ball) is often used for strength training. It is used for rehabilitation from an injury and therefore plays an important role in the field of **sports medicine**.

Moreover, it is the ideal attribute for training multiple muscle groups because of the combination of cardio and strength training.

Medicine ball training is one of the oldest forms of strength and conditioning training.

### SPECIFICATIONS

Rack measurements: 28 x 60 x 173 cm (L x W x H)

Reference code	Description	Weight	Ø	Colour	Unit
O6-O3OIO3	Mambo Max Medicine Ball   <b>Yellow</b>	1 kg	20 cm		1 pc
O6-O3OIO4	Mambo Max Medicine Ball   <b>Red</b>	1,5 kg	20 cm		1 pc
O6-O3OIO5	Mambo Max Medicine Ball   <b>Green</b>	2 kg	20 cm		1 pc
O6-O3OIO7	Mambo Max Medicine Ball   <b>Blue</b>	3 kg	22 cm		1 pc
O6-O3OIO8	Mambo Max Medicine Ball   <b>Black</b>	4 kg	22 cm		1 pc
O6-O3OIO9	Mambo Max Medicine Ball   <b>Grey</b>	5 kg	27 cm		1 pc
O6-O0O3OI	Mambo Max Medicine Ball Rack ( <b>empty</b> )				1 pc



## Kettlebell

The Mambo Max Kettlebells are an important part of any fitness or CrossFit training. Made from excellent iron surrounded with a colourful vinyl coating, they will not damage any surfaces when being stored. Starting from 2 kg up to 24 kg, we also offer you a complete set, including a premium quality Kettlebell Rack. This rack can hold up to 150 kg per shelf (total of 300 kg).

### SPECIFICATIONS

Rack measurements: 120 x 50 x 70 cm (L x W x H)



Reference code	Description	Weight	Colour	Unit
O6-O5OIO2	Mambo Max Kettlebell   <b>Sky Blue</b>	2 kg		1 pc
O6-O5OIO4	Mambo Max Kettlebell   <b>Purple</b>	4 kg		1 pc
O6-O5OIO6	Mambo Max Kettlebell   <b>Yellow</b>	6 kg		1 pc
O6-O5OIO8	Mambo Max Kettlebell   <b>Red</b>	8 kg		1 pc
O6-O5OIO10	Mambo Max Kettlebell   <b>Green</b>	10 kg		1 pc
O6-O5OIO12	Mambo Max Kettlebell   <b>Blue</b>	12 kg		1 pc
O6-O5OIO16	Mambo Max Kettlebell   <b>Black</b>	16 kg		1 pc
O6-O5OIO20	Mambo Max Kettlebell   <b>Silver</b>	20 kg		1 pc
O6-O5OIO24	Mambo Max Kettlebell   <b>Orange</b>	24 kg		1 pc
O6-O5IOOI	Mambo Max Kettlebell Rack ( <b>empty</b> )			1 pc
O6-O5IOOI	Mambo Max Kettlebell Set + Rack			Set



## Wrist & ankle weights

To quickly build stamina and muscle, simply add Mambo Max Wrist and Ankle Weights to your workout. They are filled with iron sand and they have a soft and comfortable nylon finishing. These Mambo Max Wrist and Ankle Weights will effectively shape and tone your body and you can wear them anywhere! Wear them around the office, out shopping or at the gym.



0,5 - 2 KG

Reference code	Description	Weight	Unit
O6-O2O1O2	Mambo Max Wrist & Ankle Weights	0,5 kg	Pair
O6-O2O1O3	Mambo Max Wrist & Ankle Weights	1 kg	Pair
O6-O2O1O4	Mambo Max Wrist & Ankle Weights	1,5 kg	Pair
O6-O2O1O5	Mambo Max Wrist & Ankle Weights	2 kg	Pair
O6-O2O1O7	Mambo Max Wrist & Ankle Weights	3 kg	Pair
O6-O2O1O8	Mambo Max Wrist & Ankle Weights	4 kg	Pair
O6-O2O1O9	Mambo Max Wrist & Ankle Weights	5 kg	Pair



3 - 5 KG

## New: Premium wrist & ankle weights

Complete your gym set and enhance your work out with these trendy wrist and ankle weights. With the premium Mambo Max wrist and ankle weights you can easily train your arm and leg muscles. You can build up the weights or simply challenge yourself during exercising. Available in 3 different weights, each with a particular colour:

- Light grey: 0,5 kg
- Grey: 1 kg
- Anthracite: 1,5 kg

### SPECIFICATIONS

Available in 3 different weights: 0,5 kg | 1 kg | 1,5 kg

Pair of the chosen weight comes in a colourbox

Description	Reference code	Unit
Wrist & ankle weights   0,5 kg   Light grey	O6-O6O1O1	pair
Wrist & ankle weights   1 kg   Grey	O6-O6O1O2	pair
Wrist & ankle weights   1,5 kg   Anthracite	O6-O6O1O3	pair



## Thumb lock wrist weight

The Mambo Max Thumb Lock Wrist Weight enables you to fix the weight to your thumb, ensuring a fixed position during your entire workout. Whether you're boxing, running or doing Zumba, the weights will remain in place!

Reference code	Description	Weight	Unit
O6-O2O2O2	Mambo Max Thumb Lock Wrist Weights	0,5 kg	Pair





# O7 Thermal





## Para Pro - paraffin bath

The Para Pro is a professional paraffin wax heating bath. Treatment with the Moves Para Pro Bath reduces arthritis in the joints, and it also relaxes the muscles. Additionally, the paraffin wax provides deep hydration and nourishment, quickly leaving your hands and feet soft and healthy.

### SPECIFICATIONS | PARA PRO

- Outer dimensions: 39 x 25 x 18 cm (L x W x H)
  - Inner dimensions: 32,5 x 18 x 13 cm (L x W x H)
  - Weight: 3,1 kg
  - 220/230 V: No transformer needed, which will lengthen the durability
  - Low power consumption once the paraffin has melted
  - 24-hours a day fully operational
  - Fast melting between 2-3 hours !
  - Adjustable temperature (40 - 60°C)
  - Digital thermometer
- Comes with UK or EU plug!



### SPECIFICATIONS | PARAFFIN PEARLS

- White fully refined paraffin wax
  - Congealing point: 53°C
  - Melting point: +/- 52-54°C
  - Odourless
  - Meets the purity requirements of BGA1, Ph.EUR.2 and FDA3
1. Bundesgesundheitsamt recommendation No. XXV B II
  2. European Pharmacopoeia, current edition
  3. Food and Drug Administration, § 172.886 Petroleum Wax



PARAFFIN PEARLS


### SPECIFICATIONS | PARAFFIN BLOCKS

- Green fully refined paraffin wax
- Melting point: 46°C
- Tea tree scented



HEAT HOLDING GLOVES

Includes heat holding gloves & 10 pairs of disposable gloves



PARAFFIN BLOCK

Reference code	Description	Unit
07-100101	Moves Para Pro Bath + 2,5 kg Paraffin Pearls   white   EU-plug	1 pc
07-100101 UK	Moves Para Pro Bath + 2,5 kg Paraffin Pearls   white   UK-plug	1 pc
07-100206	Moves Para Pro Bath + Set of 6 Paraffin Blocks 46°C (2,72 kg)   tea tree green   EU-plug	1 pc
07-100206 UK	Moves Para Pro Bath + Set of 6 Paraffin Blocks 46°C (2,72 kg)   tea tree green   UK-plug	1 pc
07-110102	Moves Paraffin Pearls - 2,5 kg   52-54 dab	1 pc
07-110125	Moves Paraffin Pearls - 25 kg   52-54 dab	1 pc
07-110206	Moves Paraffin Blocks 46°   Set of 6 x 453 g (2,72 kg)   tea tree green	1 pc
07-109801	Moves Insulating gloves for Para Pro Bath   blue	1 pair
07-109802	Moves Disposable Gloves for Para Pro Bath	10 pairs



**Coming soon:**

## Hydrocollator and moist heat packs

Constructed of seamless, stainless steel, a hydrocollator is a stationary or mobile thermostatically controlled liquid heating device. Once heat packs are inserted into the water inside, the hydrocollator heats them up to maximum 76°C. From here, the health professional removes the pack to wrap it in multiple towels and apply to the area being treated. The packs themselves draw in the heat to hold their temperature.

This treatment assists with pain management, addressing spasms, improving a patient's range of motion and loosening tight, contracted muscles before a massage, physical therapy or athletic training session or chiropractic adjustment.

### **MOVES MOIST HEAT PACK UNITS**

The unit heats moist heat packs up in a water bath up to 76°C. We offer 4, 6 and 12 pack unit capacity.

Each Moves moist heat pack unit is made from durable seamless, stainless steel which makes them easy to clean. All heating units are fully insulated. This insulation prevents the exterior from becoming too hot and it keeps the temperature of the heating unit inside.

### **MOIST HEAT PACKS**

The moist heat packs are easy to heat up and they retain the moist heat for up to 30 minutes. They're reusable and can be used with a cover or a towel.

### **MOIST HEAT PACK COVERS**

We also offer moist heat pack covers, to cover a hot pack unit. These are reusable and washable.



MOIST HEAT PACK IN A COVER



HYDROCOLLATOR (HEATING UNIT)



MOIST HEAT PACKS



## Hot/cold pack MD CE

The Moves Hot/Cold Packs are used for hot as well as cold therapy. This dual applicability makes them a strong alternative to ordinary heating pillows and ice packs. The versatile, reusable packs deliver therapy effectively and economically. The use of the Moves Hot/Cold Pack is unlimited and its material is non-toxic. Even below -20°C it remains soft. The Moves Hot/Cold Pack can be bent in any angle, allowing maximum coverage of joint and muscle. All Hot/Cold Packs can be heated (microwave) and can be frozen (freezer). The temperature holds for approximately 30 minutes.

### CLASSIC - GEL

The classic (and economical) transparent Hot/Cold Pack. The special gel is non-toxic and biodegradable.

✓ Sold per 10 pcs

🏠 For home use only!



### STANDARD - GEL

Our standard Hot/Cold Pack is a nylon covered gelpack. The special gel is non-toxic and biodegradable.

✓ Sold per 10 pcs



### SOFT TOUCH - CLAY

The soft touch provides a user-friendly and frost-free surface that allows usage in **direct contact with the skin**. Filled with clay, a 100% natural product.


✓ Sold per 10 pcs



### PROTECTION SLEEVE

100% polyester Protection Sleeve for Hot/Cold Packs. Machine washable. Available in 3 sizes (S | M | L).

✓ Sold per 10 pcs



### TRI-SECTIONAL HOT/COLD PACK

The soft touch tri-sectional Hot/Cold Pack treats multiple zones at once (neck and shoulders).

✓ Sold individually & per 10 pcs



### XXL HOT/COLD PACK

Ideal for pain relief on shoulders and back. This Hot/Cold Pack has a nylon cover.

✓ Sold individually



Dimension		Moves Hot/Cold Packs				Unit
		Classic (Gel)	Standard (Gel)	Soft Touch (Clay)	Protection Sleeve	
Small	15 x 25 cm	07-010102	07-010202	07-010302	07-020102	10 pcs
Medium	20 x 30 cm	07-010103	07-010203	07-010303	07-020103	10 pcs
Large	25 x 35 cm	07-010104	07-010204	07-010304	07-020104	10 pcs
XXL	33 x 47 cm	n/a	07-010206	n/a	n/a	1 pc
Tri-Sectional	40 x 20 cm	n/a	07-010211	07-010311	n/a	10 pcs
Tri-Sectional	40 x 20 cm	n/a	n/a	07-010312	n/a	1 pc



## PearlPack

Hot and cold therapy can be used to relieve pain and inflammations associated with muscle strain and sprain. This dual applicability makes them a strong alternative to ordinary heating pillows and ice packs. The small pearls absorb heat or cold, which they can retain for 20 minutes. During this time, these packs will relieve a part of your pain.

Not only do they look good, the products are very convenient in use. PearlPacks won't feel wet, nor drip after freezing or heating them, avoiding a mess. You can also use these in **direct contact with the skin**. All materials used are non-toxic and odourless.

Besides relieving pain from an injury or muscle strain, hot/cold packs are often used to reduce swelling after an operation. PearlPacks truly are multifunctional products.

### PEARLPACK | SPORTS PACK & STRAP

The PearlPack Sports Pack comes with a handy velcro strap.




### PEARLPACK | BACK WRAP

This large pack divided in 3 segments is equipped with a velcro strap. It fits perfectly around the user's waist for an optimal pain relief on the back. The pearls can be pushed from one segment to the other if one area of the back needs less or more thermal treatment.



### PEARLPACK | EYE MASK

This mask adapts to the user's face and has a fabric strap. Hot treatment can be applied when the eyes feel dry or if there's an infection of the eyelid. Cold treatment will help when suffering from headaches, swollen eyes or sinusitis.



### PEARLPACK | KNEE/ELBOW

This offers a solution for hot/cold therapy on the knee or elbow joints. The large opening ensures that your knee or elbow remains free, so you can continue to move freely, while the Pearl technology sends direct relief to the pain area. The flexible velcro strap keeps the pack firmly in place.






Reference code	Description	Dimension	Unit
07-O40102	Moves PearlPack   Sports Pack+ strap	19 x 12,5 cm	1 pc
07-O40201	Moves PearlPack   Back Wrap	43 x 17,5 cm	1 pc
07-O41101	Moves PearlPack   Eye Mask	22,5 x 7,5 cm	1 pc
07-O40301	Moves PearlPack   Knee/Elbow	37 x 25,6 cm	1 pc



**Roll-on pack** MD CE

Moves Roll-On Packs combine hot/cold treatment with practicality. The tubular shape and flexible material allow an easy roll-on application, even after being heated or frozen. The various sizes will fit perfectly over your finger, wrist, ankle or knee. This design ensures that the whole area will benefit from the treatment.

**SPECIFICATIONS**

-  Available in 5 sizes
-  Reusable, non-toxic and latex-free
-  Each Roll-On Pack includes a protective sleeve



Reference code	Description	Inside diameter Ø	Circumference C	Uses	Unit
O7-O5O1O1	Moves Roll-On Pack   Finger	2,5 cm	One size fits all	Finger	1 pc
O7-O5O1O2	Moves Roll-On Pack   Small	10,0 cm	Fits up to 25 cm	Hand, wrist, arm, elbow	1 pc
O7-O5O1O3	Moves Roll-On Pack   Medium	12,5 cm	Fits 25 - 38 cm	Hand, wrist; arm, elbow, ankle	1 pc
O7-O5O1O4	Moves Roll-On Pack   Large	15,0 cm	Fits 38 - 53 cm	Elbow, ankle, foot, calf, knee	1 pc
O7-O5O1O5	Moves Roll-On Pack   Extra Large	17,5 cm	Fits 53 - 58 cm	Knee, thigh	1 pc




## Instant heat pack

Have you ever needed a heat pack, but there was no kitchen or microwave in sight? Moves has the perfect solution for you! Your heat pack on-the-go. Simply flex the metal disk inside our Moves Instant Heat Pack and **within seconds** it will warm up to around a comfy **54°C**. Furthermore, it is **reusable** many times! Just boil the pack for 10 to 15 minutes and it will return to its original fluid state.

The Moves Instant Heat Pack comes in **two different sizes**:


### 1. POCKET SIZE PACK

Perfect for outdoor activities! Easy to use while camping, going to sporting events or whenever you go outside in colder weather. It is small, lightweight and fits neatly in your pocket.



### 2. SHOULDER & NECK PACK

This 3-piece instant heat pack is designed to fit comfortably around your neck and shoulders. The instant heat pack loosens your muscles and relieves pressure and soreness.




Reference code	Description	Dimension	Unit
O7-O6OIOI	Moves Instant Heat Pack   Pocket   Set of 3	12,7 x 7,6 cm	1 pc
O7-O6OIII	Moves Instant Heat Pack   Shoulder & Neck	41 x 25 cm	1 pc





# O8 Evaluation






## New: MicroFET™2 muscle test dynamometer



The wireless microFET™2 Digital Handheld Dynamometer is an accurate, portable Force Evaluation and Testing (FET) device. It is designed for taking objective, reliable and quantifiable **muscle testing measurements**. The microFET™2 aids in diagnosis, prognosis and treatment of neuromuscular disorders.

- Measurement range 0-136 kg (0-300 lb) force
- Easy-to-read LCD Displays shows peak force and elapsed time
- Selectable units of measure: pounds (lbs.), Newtons (N), or kilogram-force (kgf)
- Accuracy within 1% of the reading
- Two threshold settings for muscle testing: Low 0,4 to 136 kg (0,8 to 300 lb) in 0,05 kg (0,1 lb) increments and High threshold 1,4 to 136 kg (3,0 to 300 lb) in 0,05 kg (0,1 lb) increments
- Stores up to 30 tests
- Uses rechargeable lithium ion rechargeable battery
- Self-activating "sleep" mode after 3 minutes to extend battery life



Description	Reference code	Unit
MicroFET™2 Muscle Test Dynamometer	12-O381W	1 pc

## Fingertip pulse oximeter



A Pulse Oximeter is a very important and common device to check the patient's **blood-oxygen saturation (SpO<sub>2</sub>) level** and **pulse rate**. As a compact, reliable and durable monitoring device, the Fingertip Pulse Oximeter greatly enhances patient care.

### JUNIOR OXIMETER

Designed for kids from 3 - 12 years old.

- Rechargeable via USB-cable
- Dual colour OLED-screen displays SpO<sub>2</sub>, PR, pulse bar and waveform
- 4-direction & 6-mode display for convenient reading
- Set-up alarm range of SpO<sub>2</sub> and pulse rate
- Menu-function setting (beep sounds, etc.)
- Lithium battery; automatically power off
- Gravity function, automatic rotating reading option



### ADULT OXIMETER

Features:

- Dual colour OLED-screen displays SpO<sub>2</sub>, PR, pulse bar and waveform
- SpO<sub>2</sub>; measurement range between 70%-99% accuracy; ± 2% on the stage of 80%-99%
- PR; measurement range between 30BPM-240BPM accuracy; ± 1BPM or ±1% (the larger one)
- Low power indicator
- 6 display modes with adjustable brightness
- Automatically power off
- AAA alkaline batteries



Reference code	Description	Unit
O8-O7O2O1	Fingertip Pulse Oximeter   Adult   OLED   2 x AAA	1 pc
O8-O7O2O2	Fingertip Pulse Oximeter   Junior 3-12y old   OLED   Rechargeable	1 pc



## Blood pressure monitor MD CE

The blood pressure monitor is used for measuring **blood pressure** and **pulse rate**. It is suitable for people older than 12 years. It features:

- Large LCD display: 72 x 43,5 mm
- 2 users, 99 memories per user
- Irregular heart beating sign
- Indicating blood pressure level in 4 colours
- Average of last 3 records
- Low power detection
- Over pressure protection
- Automatic switch-off



Reference code	Description	Unit
O8-O7O5O1	Blood Pressure Monitor + Cuff	1 pc

## Infrared thermometer MD CE

### EAR & FOREHEAD THERMOMETER

The infrared thermometer is used for measuring body temperature by two-in-one functions: Ear and Forehead temperature. It is suitable for both medical and home use.

- Fast and easy temperature measurement by advanced infrared technology
- Range between 34,0°C-43,0°C accuracy; <35,0°C and >42,0°C = ± 0,3°C
- Memory mode
- Suitable for both infants and adults
- Easy operating and cleaning
- Automatic power off after 60 seconds
- 2 AAA alkaline batteries
- Tri-colour backlit:
  - <37,4°C
  - 37,4°C ~ 38,0°C
  - >38,0°C



### SHOTGUN THERMOMETER

Measures your body temperature without direct contact with the skin. Can measure between 5-15 cm from the skin.

- Detailed 34 memories with time/data
- Response time: 1 second
- Effective distance: 5-15 cm
- User guidance by infrared beam
- Fever alarm



Reference code	Description	Unit
O8-O7O3O1	Infrared Thermometer   Ear & Forehead auto-sensing   2 x AAA	1 pc
O8-O7O3O3	Infrared Thermometer   Shotgun   Forehead auto-sensing   2 x AAA	1 pc



## Pedi-scoliometer



Scoliometers are inclinometers designed to **measure trunk asymmetry or Angle of Truncal Rotation (ATR)** used in scoliosis screening programmes. They are small, non-invasive devices which are placed over the spine while a person is in the forward bending position. The scoliometer reading indicates the degree of trunk asymmetry (rotation).

Reference code	Description	Unit
08-060102	Pedi-Scoliometer	1 pc



## moves

### Neurological reflex instruments



#### 1. BUCK HAMMER

Two heads for cutaneous sensation evaluation. Comes with a needle and a brush.

#### 2. BABINSKI HAMMER

Pointed bottom can be used as gross aesthesiometer.

#### 3. TAYLOR HAMMER

Used to elicit responses during an examination.

#### 4. DEJERINE HAMMER

Used for testing deep tendon reflexes in a diagnostic evaluation. Double-sided head composed of two rounded rubber surfaces for increased comfort.


#### 5. TROEMNER HAMMER

Shaped like a two-headed mallet. The larger mallet is used to elicit tendon stretch reflexes, and the smaller mallet is used to elicit percussion myotonia.


#### 6. WARTENBERG PINWHEEL

Used to elicit cutaneous sensory and pain perception responses.


Reference code	Description	Unit
08-050101	Moves Buck Hammer	1 pc
08-050102	Moves Babinski Hammer	1 pc
08-050103	Moves Taylor Hammer	1 pc
08-050104	Moves Dejerine Hammer	1 pc
08-050105	Moves Troemner Hammer	1 pc
08-050201	Moves Wartenberg Pinwheel	1 pc



1. BUCK HAMMER




2. BABINSKI HAMMER



3. TAYLOR HAMMER



4. DEJERINE HAMMER



5. TROEMNER HAMMER



6. WARTENBERG PINWHEEL



# <sup>®</sup>**SAEHAN**

Grip



## **About the Saehan® brand**

Saehan® provides a wide range of measurement products. For example dynamometers, goniometers, body fat calipers and so on.

All Saehan® products are CE and FDA approved.

Discover the product assortment [here!](#)



## Hydraulic hand dynamometer MD CE

The Hydraulic hand dynamometer gives accurate grip strength readings without the patient being able to “feel” the handle move. The adjustable handle has 5 positions to adapt to any hand size. The indicator remains at the patient’s maximum reading until it is reset. This product has a scale in both kilograms and pounds.

### SPECIFICATIONS

- ↔ Range: 0 - 90 kg (200 lb)
- ✓ Comes in a quality hard carrying case
- ✓ Warranty: 2 years

Reference code	Description	Unit
O8-O1O1O1	Saeahan Hydraulic Hand Dynamometer	1 pc



## Hydraulic 3-piece hand evaluation kit MD CE

A portable hard carrying case containing all instruments necessary to evaluate the strength and range-of-motion of the hand and fingers.

### SPECIFICATIONS

- ✓ Includes: hydraulic hand dynamometer
- ✓ Includes: hydraulic pinch gauge
- ✓ Includes: stainless steel finger goniometer 15 cm - 180° per 5°
- ✓ Warranty: 2 years

Reference code	Description	Unit
O8-O1O113	Saeahan Hydraulic Hand Evaluation Kit   3 pieces	Set



## Hydraulic 7-piece hand evaluation kit MD CE

This 7-piece Hand Evaluation Kit is a complete set of instruments to evaluate the strength and range-of-motion of the hand and fingers.

### SPECIFICATIONS

- ✓ Includes: hydraulic hand dynamometer
- ✓ Includes: hydraulic pinch gauge
- ✓ Includes: stainless steel finger goniometer 15 cm - 180° per 5°
- ✓ Includes: 2-point discriminator with 3<sup>rd</sup> point
- ✓ Includes: Wartenberg pinwheel
- ✓ Includes: finger gauge
- ✓ Includes: functional finger motion gauge
- ✓ Warranty: 2 years

Reference code	Description	Unit
O8-O1O117	Saeahan Hydraulic Hand Evaluation Kit   7 pieces	Set





## Digital hand dynamometer with data transfer software - DHD-1 CE

Designed for reliable grip-force evaluations. This hand dynamometer has an easy-to-read LCD display. It features an electronic zero calibration system and up to 40 test data can be stored in memory. It is powered by a rechargeable Li-Ion polymer battery and data will not be lost when changing the battery. The Data Transfer Software can transfer the test results into .txt-files, which can be easily accessed using Notepad or MS-Word for editing and printing.

### SPECIFICATIONS

- ↔ Range: 0 - 90 kg (200 lb)
- ✓ Comes in a quality hard carrying case
- ✓ Includes: digital hand dynamometer
- ✓ Includes: DTS-software
- ✓ Includes: USB cable
- ✓ Warranty: 1 year



Reference code	Description	Unit
08-010201	Saehan Digital Hand Dynamometer - DHD-1	1 pc


## Digital hand dynamometer with G-STAR™ software - DHD-3 CE

The DHD-3 Digital Hand Dynamometer (SH1003) provides the same features as the DHD-1 when used as a stand-alone device, but adds automated Rapid Exchange Grip Testing. Using Auto Mode, up to 10 strength tests for each hand, as well as the time allowed to switch hands, can be configured to automate the testing procedure.

The Grip Strength Testing and Research (G-STAR™) is an easy to use software tool which automates grip strength testing and endurance testing in both clinical and research environments. All patient information and their test results are stored in an Access database.

### SPECIFICATIONS

- ↔ Range: 0 - 90 kg (200 lb)
- ✓ Comes in a quality hard carrying case
- ✓ Minimum requirements: Windows XP and USB port
- ✓ Warranty: 1 year



Reference code	Description	Unit
08-010202	Saehan Digital Hand Dynamometer - DHD-3 (including G-STAR™ software)	1 pc

## Squeeze dynamometer MD CE

The Squeeze (bulb) Dynamometer is an inexpensive device to reliably measure hand and finger strength. The dynamometer has an indicator dial that remains at the maximum reading until reset. Unit is calibrated in pounds and kilograms.

### SPECIFICATIONS

- ↔ Range: 0 - 70 kg (150 lb)
- ✓ Comes in a quality soft pouch case
- ✓ Warranty: 1 year



Reference code	Description	Unit
08-010301	Saehan Squeeze Dynamometer   70 kg	1 pc



## Smedley hand dynamometer

MD CE

The Smedley spring-type Dynamometer has been successfully used for many years to measure grip strength. The ergonomic handle is adjustable to fit all hand sizes. The unit is calibrated in kilograms and pounds.

### SPECIFICATIONS

- ↔ Range: 0 - 100 kg (220 lb)
- ✓ Comes in a quality soft pouch case
- ✓ Warranty: 1 year



Reference code	Description	Unit
08-010401	Saehan Smedley Hand Dynamometer   Ergonomic handle	1 pc

## Back-leg-chest dynamometer


MD CE

This Dynamometer measures the strength of your back, leg and chest. The base provides secure footing. The chain length is adjustable to accommodate for height differences or to vary the point of force application. It shows both pounds and kilograms. The pointer remains at maximum until it is reset.

### SPECIFICATIONS

- ↔ Range: 0 - 300 kg (660 lb)
- ✓ Comes fully assembled + 1,5 m chain
- ✓ Warranty: 1 year

Reference code	Description	Unit
08-010601	Saehan Back-Leg-Chest Dynamometer	1 pc



## Pinch gauge

MD CE

### MECHANICAL PINCH GAUGE

Uses for hand-strength testing by measuring tip, key and palmar pinch force. Made of durable painted stainless steel. Mechanical properties will not degrade over time. It can serve as a diagnostic tool or for monitoring the progress of rehabilitation



### SPECIFICATIONS

- ↔ Range: 0 - 20 kg (45 lb)
- ✓ Warranty: 1 year

Reference code	Description	Unit
08-020102	Saehan Mechanical Pinch Gauge	1 pc

### HYDRAULIC PINCH GAUGE

Uses the improved hydraulic system of the hand dynamometer to assure convenience, product reliability, measurement accuracy and repeatability. This results in a more accurate test result for all pinch tests.



### SPECIFICATIONS

- ↔ Range: 0 - 20 kg (45 lb)
- ✓ Warranty: 2 years

Reference code	Description	Unit
08-020101	Saehan Hydraulic Pinch Gauge	1 pc



## Plastic goniometer




Clear plastic permits observation of joint's axis of motion, and its range of motion. The 360° head had three scales calibrated to be used with ISOM (International Standard Of Measurements, STFR) system.




15 CM GONIOMETER




20 CM GONIOMETER



30 CM GONIOMETER




RULONG GONIOMETER



POCKET GONIOMETER

Reference code	Description	Length	Unit
08-030101	Saehan Plastic Goniometer   0° to 360° per 1°	15 cm	1 pc
08-030102	Saehan Plastic Goniometer   0° to 360° per 1°	20 cm	1 pc
08-030103	Saehan Plastic Goniometer   0° to 360° per 1°	30 cm	1 pc
08-030111	Saehan Rulong Goniometer   0° to 360° per 2°	18 cm	1 pc
08-030112	Saehan Pocket Goniometer   180° x 2 per 5°	15 cm	1 pc


## Stainless steel goniometer




08-030201



08-030202



08-030203  
ARM TENSION IS CONTROLLED BY THUMB KNOB



08-030204  
EXTRA LONG ARMS FOR MEASURING ALL JOINTS


Reference code	Description	Length	Scale	Increments	Unit
08-030201	Saehan Stainless Steel Robinson Pocket Goniometer	15 cm	180°	5°	1 pc
08-030202	Saehan Stainless Steel Goniometer	20 cm	180° x 2	1°	1 pc
08-030203	Saehan Stainless Steel Goniometer	35 cm	180° x 2	1°	1 pc
08-030204	Saehan Stainless Steel Goniometer	35 cm	360°	1°	1 pc



## Finger goniometer

MD CE

Finger and Flexion/Hyper-Extension Goniometers measure range-of-motion of finger joints. The head has two opposing 180° scales marked in 5° increments. A linear scale in centimeters is located on the arm.




Reference code	Description	Unit
O8-O3O3O1	Saeahan Stainless Steel Finger Goniometer   9 cm   180° per 5°	1 pc
O8-O3O3O2	Saeahan Stainless Steel Finger Goniometer   15 cm   180° per 5°	1 pc

## Goniometer kit

MD CE

Excellent soft pouch case!



### PLASTIC SET

- Plastic Goniometer 15 cm | 0° to 360° per 1°
- Plastic Goniometer 20 cm | 0° to 360° per 1°
- Plastic Goniometer 30 cm | 0° to 360° per 1°
- Rulong Goniometer 18 cm | 0° to 360° per 2°
- Plastic Goniometer 15 cm | 180° x 2
- Plastic Flexion-Hyperextension Goniometer



### STAINLESS STEEL SET


- Robinson Pocket Goniometer 15 cm | 180°
- Stainless Steel Goniometer 18 cm | 180°
- Stainless Steel Goniometer 35 cm | 180° x 2
- Stainless Steel Goniometer 35 cm | 360°
- Stainless Steel Finger Goniometer 15 cm
- X-Ray Goniometer 20 cm | 180°

Reference code	Description	Unit
O8-O3O121	Saeahan Goniometer Set   Plastic   6 pcs	Set
O8-O3O211	Saeahan Goniometer Set   Stainless Steel   6 pcs	Set

## Tuning fork

These aluminium Tuning Forks are developed to measure sensitivity to vibration which can be an indicator of nerve damage. It is also used for foot exams by diabetic people to check for signs of neuropathy.

Reference code	Description	Unit
O8-O5O3O1	Tuning Fork   128cps	1 pc
O8-O5O3O2	Tuning Fork   256cps	1 pc





## Fat caliper MD CE

Professional light-weight aluminium caliper. The measurement arms have floating tips to assure accurate and reproducible results. The 65 mm scale is large and easy to read. A corrective scale for sex and age is included in a multilingual manual

### SPECIFICATIONS

- ✓ Comes in a quality soft pouch carrying case

Reference code	Description	Unit
08-040101	SaeHan Skinfold Fat Caliper Aluminium - Green	1 pc







### **About the Baseline® brand**

Baseline® measurement instruments have been used by orthopedists, therapists, FCE administrators and professors since the 1980s. Many of you who have received therapy after an injury have probably had your progress tested throughout your rehabilitation with a Baseline® measurement instrument.

The product line includes instruments to measure strength, range-of-motion (ROM), sensory perception, body measurements and more.



## It's time to upgrade to digital with the BIMS™ digital 5-position grip and pinch dynamometers!

### WHY SWITCH FROM A HYDRAULIC DYNAMOMETER TO DIGITAL?

- Force is measured using load-cell digital technology that gives more accurate, repeatable, and consistent measurements
- No hydraulic oil is used to measure force, eliminating the risk of oil leakage and decalibration
- Resetting of max results is digital; standard hydraulic reset knobs are prone to breaking
- LCD screen makes results easy to read
- Saves time by calculating mean, SD, COV, and other statistics on the device

### WHY SHOULD I CHOOSE BIMS™ DIGITAL DYNAMOMETERS?

- Extended testing capacity (135 kg grip (300 lb) / 45 kg (100 lb) pinch) is designed for the strongest and weakest patients
- Measures force consistently with  $\pm 1\%$  accuracy
- Color LCD display graphically shows results
- Numerous testing protocols for advanced testing
- 5-position hand grip accommodates different hand sizes and tests
- Results saved on device with memory for over 1000 patients





## BIMS™ digital (Load-Cell) 5-position dynamometers

- Adjustable 5-position handle fits all size hands
- Accurate and repeatable grip and pinch strength measurements
- Grip dynamometers have an extended testing capacity that measures grip strength up to 135 kg (300 lbs)
- Pinch dynamometers measure pinch strength up to 45 kg (100 lbs)
- Purchase the dynamometer that meets your practice needs: **Functional, Clinic** or **Deluxe**
- Subject/Test information can be entered at time of test: for example, dominant hand, handle position, exertion/rest times, units (lb/kg)
- **1 year warranty**

### TEST PROTOCOLS

- **Live Test:**  
Shows real time and maximum value. Test result is not stored
- **Quick Test:**  
Shows real time and maximum value. Maximum is stored
- **GST (General Strength Test):**  
3 tests are performed at the same handle position. Peak, mean, SD, and COV are stored
- **RET (Rapid Exchange Test):**  
6 tests (alternating hands) are performed at the same handle position
- **MMVE (Modified Maximum Voluntary Effort Test):**  
5 tests are performed; 1 at each handle position
- **MVE (Maximum Voluntary Effort Test):**  
15 tests are performed per hand; 3 at each handle position
- **Fatigue (Work Test):**  
1 extended time test is performed. Results can be used to compare strength exerted (total work performed) over 2 or more user defined time periods



3-PIECE SETS INCLUDE GRIP DYNAMOMETER, PINCH DYNAMOMETER, PLASTIC FINGER GONIOMETER (12-1014), CARRY CASE, INSTRUCTIONS AND NORMS

### FUNCTIONAL MODEL MD CE

This model does everything a standard hydraulic dynamometer can do with greater accuracy and repeatability. Use to view **live strength readings** with the **Max Force** test.

### CLINIC MODEL MD CE

This model is ideal for clinic practitioners to perform timed tests which are **calculated and saved on the device**. It includes the **Max Force**, **Quick**, and **GST** tests.

### DELUXE MODEL MD CE

This model expands upon the Clinic model and includes **more advanced timed tests** that are ideal for research applications, Workers' Compensation, FCE (Functional Capacity Evaluation) and day-to-day testing. It includes **Max Force**, **Quick**, **GST**, **RET**, **MMVE**, **MVE**, and **Fatigue** tests.

Reference code	Description	Unit
12-0072	Grip (135 kg / 300 lb)	1 pc
12-0082	Pinch (45 kg / 100 lb)	1 pc
12-0092	3-piece hand set	set

Reference code	Description	Unit
12-0070	Grip (135 kg / 300 lb)	1 pc
12-0080	Pinch (45 kg / 100 lb)	1 pc
12-0090	3-piece hand set	set

Reference code	Description	Unit
12-0071	Grip (135 kg / 300 lb)	1 pc
12-0081	Pinch (45 kg / 100 lb)	1 pc
12-0091	3-piece hand set	set



## 5-position hydraulic hand grip dynamometer

MD CE

- This product gives accurate grip strength readings without the subject being able to "feel" the handle move
- Readings are repeatable
- Maximum reading remains until reset
- Readings in lb and kg
- CE certified; manufactured and serviced in the USA since 1978
- Recalibration available once/year based on warranty
- Five position handle and body contours assure results consistent with published Baseline® and Jamar® studies
- Internationally accepted design assures reliability, accuracy and repeatability
- Lightweight units include protective carry case, instructions for use and norms
- **90 kg (200 lb) capacity, 6,4 cm (2,5") diameter dial head**

### CHOOSE YOUR WARRANTY:



LITE® - 1 YEAR WARRANTY



STANDARD - 2 YEARS WARRANTY



HD® - 5 YEARS WARRANTY

Reference code	Description	Colour	Warranty	Unit
12-O241	LITE®   Red finish	Red	LITE®   1 year warranty	1 pc
12-O240	Standard   Blue finish	Blue	Standard   2 year warranty	1 pc
12-O221	HD®   Black finish	Black	HD®   5 year warranty	1 pc

## Pneumatic squeeze dynamometers

MD CE

- The inexpensive way to measure grip and pinch strength
- Available with and without a maximum force indicator (reset) that remains at the maximum reading until reset
- Fluid-filled bulb; calibrated in PSI
- Smaller bulb 2,5 cm (1") can be used for children and pinching
- Often used in pediatric rehabilitation



Reference code	Description	Unit
12-O291	30 PSI   1 year warranty	1 pc
12-O293	15 PSI   1 year warranty	1 pc

5 CM (2") BULB WITH RESET

Reference code	Description	Unit
12-O297	3-bulb set   2 year warranty	set

3 BULB SET FOR GRIP AND PINCH MEASUREMENTS, INCLUDES 2,5 CM (1"), 3,8 CM (1,5") AND 5 CM (2") BULBS WITH RESET AND CASE. SIMILAR TO THE MARTIN VIGOROMETER



## Smedley spring dynamometer

MD CE

- Adjustable handle to fit all hand sizes
- Readings in both kg and lb
- Includes carrying case
- Ideal cost-efficient solution to measure grip strength
- **2 years warranty**



Reference code	Description	Unit
12-O281	Adult (100 kg / 220 lb)	1 pc

## Electronic Smedley dynamometer

MD CE

- 5-preset grip sizes
- 90 kg (200 lb) capacity
- Saves/stores results for up to 19 users for fast retrieval
- **2 years warranty**



Reference code	Description	Unit
12-O286	Digital Smedley (90 kg / 200 lb)	1 pc



## Mechanical pinch gauge

MD CE

- Measure tip, key, and palmar pinch strength in kg and lb
- Results are consistent with published Mathiowetz studies
- Indicator remains at the maximum reading until reset
- Includes protective case
- **2 years warranty**



Reference code	Description	Colour	Purpose	Unit
12-O200	13,6 kg (30 lb), blue	Blue	For orthopedic and sports medicine	1 pc
12-O201	27 kg (60 lb), red	Red	For orthopedic and sports medicine	1 pc
12-O203	4,5 kg (10 lb), silver	Silver	For weak or damaged hands	1 pc
12-O202	0,9 kg (2 lb), gold	Gold	For weak or damaged hands	1 pc





## Hydraulic pinch gauge MD CE

- Measure tip, key, and palmar pinch strength in kg and lb
- Includes protective case
- The hydraulic system guarantees a more accurate test result for all pinch tests

### CHOOSE YOUR WARRANTY:

Reference code	Description	Colour	Warranty	Unit
12-O226	23 kg (50 lb) LITE®	Red	LITE®   1 year warranty	1 pc
12-O235	23 kg (50 lb) Standard	Blue	Standard   2 year warranty	1 pc
12-O222	23 kg (50 lb) HD®	Black	HD®   5 year warranty	1 pc




## 5-position hydraulic pinch gauge MD CE

Unlike standard pinch gauges, the 5-position pinch gauge has an easy-to-adjust paddle that accommodates to any hand size. This assures **accurate and repeatable pinch measurements**. 5 pinch positions permit strength test protocols (MVE and MMVE) to be used with tip, palmar and key tests.

- Users can't see results during testing
- 23 kg (50 lb) capacity
- Reads in both kg and lb
- **2 years warranty**

Reference code	Description	Unit
12-O480	5-position hydraulic pinch gauge	1 pc



## Electronic pinch gauge MD CE

- LCD display shows grip force in kg or lb
- Perform tip, key, and palmar tests
- Stores 99 readings
- 27 kg (60 lb) capacity
- 3V battery included
- **2 years warranty**

Reference code	Description	Unit
12-O475	Baseline® electronic pinch gauge	1 pc





## 7-piece hand sets



- Protective carrying case stores instruments
- 3 different sets based on the **hydraulic hand grip dynamometer of your choice**.



### 1. LITE® SET (12-O102)

This set contains the following components:

- **Hydraulic hand grip dynamometer LiTE® | Red finish | 1 year warranty (12-O241)**
- Mechanical pinch gauge | 13,6 kg (30 lb) | Blue finish (12-O200)
- Stainless steel goniometer
- 2-point discriminator with 3<sup>rd</sup> point
- Wartenberg pinwheel
- Finger circumference gauge
- Functional finger motion gauge

Reference code	Description: Grip / Pinch	Warranty	Unit
12-O102	90 kg (200 lb) (12-O241) / 13 kg (30 lb) (12-O200) LiTE®	LiTE®   1 year warranty	Set



SET 1

### 2. STANDARD SET (12-O100)

This set contains the following components:

- **Hydraulic hand grip dynamometer Standard | Blue finish | 2 year warranty (12-O240)**
- Mechanical pinch gauge | 13,6 kg (30 lb) | Blue finish (12-O200)
- Stainless steel goniometer
- 2-point discriminator with 3<sup>rd</sup> point
- Wartenberg pinwheel
- Finger circumference gauge
- Functional finger motion gauge

Reference code	Description: Grip / Pinch	Warranty	Unit
12-O100	90 kg (200 lb) (12-O240) / 13 kg (30 lb) (12-O200) Standard	Standard   2 year warranty	Set



SET 2

### 3. HD® SET (12-O117)

This set contains the following components:

- **Hydraulic hand grip dynamometer HD® | Black finish | 5 year warranty (12-O221)**
- Mechanical pinch gauge | 13,6 kg (30 lb) | Blue finish (12-O200)
- Stainless steel goniometer
- 2-point discriminator with 3<sup>rd</sup> point
- Wartenberg pinwheel
- Finger circumference gauge
- Functional finger motion gauge

Reference code	Description: Grip / Pinch	Warranty	Unit
12-O117	90 kg (200 lb) (12-O221) / 13 kg (30 lb) (12-O200) HD®	HD®   5 year warranty	Set



SET 3



## 3-piece hand sets



- Protective carrying case stores instruments
- 3 different sets based on the **warranty of your choice**:



### 1. LITE<sup>®</sup> SET (12-O101)

This set contains the following components:

- **Hydraulic hand grip dynamometer LiTE<sup>®</sup> | Red finish | 1 year warranty (12-O241)**
- Hydraulic pinch gauge | 23 kg (50 lb) | LiTE<sup>®</sup> | Red finish | 1 year warranty (12-O226)
- 15 cm (6") stainless steel goniometer

Reference code	Description Grip / Pinch	Warranty	Unit
12-O101	90 kg (200 lb) (12-O241) / 23 kg (50 lb) (12-O226) LiTE <sup>®</sup>	LiTE <sup>®</sup>   1 year warranty	Set



SET 1

### 2. STANDARD SET (12-O103)

This set contains the following components:

- **Hydraulic hand grip dynamometer Standard | Blue finish | 2 year warranty (12-O240)**
- Hydraulic pinch gauge | 23 kg (50 lb) | Standard | Blue finish | 2 year warranty (12-O235)
- 15 cm (6") stainless steel goniometer

Reference code	Description Grip / Pinch	Warranty	Unit
12-O103	90 kg (200 lb) (12-O240) / 23 kg (50 lb) (12-O235) Standard	Standard   2 year warranty	Set



SET 2

### 3. HD<sup>®</sup> SET (12-O116)

This set contains the following components:

- **Hydraulic hand grip dynamometer HD<sup>®</sup> | Black finish | 5 year warranty (12-O221)**
- Hydraulic pinch gauge | 23 kg (50 lb) | HD<sup>®</sup> | Black finish | 5 year warranty (12-O222)
- 15 cm (6") stainless steel goniometer

Reference code	Description Grip / Pinch	Warranty	Unit
12-O116	90 kg (200 lb) (12-O221) / 23 kg (50 lb) (12-O222) HD <sup>®</sup>	HD <sup>®</sup>   5 year warranty	Set



SET 3



## Push-pull dynamometer


MD CE

- Lightweight, small, and portable; easy to grasp with or without handles
- Use for muscle strength management, job task analysis, or functional capacity evaluation
- Comes with 3 push pads and 2 pull hooks
- **2 years warranty**

Reference code	Description	Unit
12-O393	Hydraulic - analog (45 kg / 100 lb)	1 pc
12-O389	Dual Grip for push-pull dynamometer	1 pc



HYDRAULIC




DUAL GRIP

## Back-leg-chest dynamometer

MD CE

- Measures strength of back, leg, and chest muscles
- Base provides secure footing
- Adjustable chain to accommodate for height differences or to vary the point of force application
- Comes fully assembled with 1,5 m (5-foot) chain
- Oversize base 60 x 60 cm (24" x 24"); regular base 45 x 45 cm (18" x 18")
- Pointer remains at maximum until reset
- Reads kg and lb

Reference code	Description	Unit
12-O400	Adult with standard base, 300 kg (660 lb)	1 pc




## Finger (digit) goniometers

MD CE

- Measures 110° flexion through 40° hyperextension in 2° increments on plastic model (5° increments on metal model)
- Arm has 10 cm (4") linear scale
- Use one finger to move the measurement head through flexion, extension, and hyperextension
- **Readings can be taken using only one hand;** freeing the therapist's other hand to stabilize the patient during measurement or to record the measurement. Standard finger goniometers require two hands to take a finger ROM measurement

Reference code	Description	Unit
12-1014	Plastic digit goniometer	1 pc
12-1016	Stainless steel digit goniometer	1 pc





## Plastic finger goniometer

MD CE

- Scale reads in 1 degree increments; linear scale reads in cm and inches
- Measures flexion to hyper-extension

Reference code	Description	Unit
12-1012	Flexion to Hyper-Extension Finger Goniometer	1 pc
(old reference: O8-O3O6O1)		




## 180° 20 cm goniometer

MD CE

- Two opposing scales marked in 1° increments
- Tension between the arms is controlled by the thumb knob
- Stainless steel

Reference code	Description	Unit
12-1040	20 cm (8"), 180°	1 pc



## 180° 35 cm goniometer

MD CE

- Two opposing scales marked in 1° increments
- Extra-long arms for measuring all joints
- Stainless steel

Reference code	Description	Unit
12-1041	35 cm (14"), 180°	1 pc




## 360° 35 cm goniometer

MD CE

- Two opposing 180° scales and one 360° scale, all marked in 1° increments
- Arm tension is controlled by the thumb knob
- Stainless steel

Reference code	Description	Unit
12-1050	35 cm (14"), 360°	1 pc





## 180° finger (small joint) goniometer

MD CE

- Stainless steel goniometer measures ROM of finger joints (metacarpophalangeal and interphalangeal) and other small joints including toes
- Head has opposing 150° scales in 5° increments
- Linear scale in cm and inches

Reference code	Description	Unit
12-1010	Standard finger, 15 cm (6")	1 pc
12-1015	Small finger, 9 cm (3,5")	1 pc



SMALL FINGER



STANDARD FINGER

## 6-piece goniometer set

MD CE

- 6 goniometers in a padded carrying case
- 15 cm (6") 180° Robinson pocket
- 35 cm (14") 180° stainless steel
- 35 cm (14") 360° stainless steel
- 15 cm (6") stainless steel finger
- 20 cm (8") 180° x-ray
- 20 cm (8") 180° stainless steel


Reference code	Description	Unit
12-1043	Goniometer set	Set




## XTender<sup>®</sup> extendable goniometer

MD CE

- Arms extend to measure large and small joints
- Extendable arms allow the goniometer to lie over anatomical bony landmarks
- Head is the same as a traditional goniometer for easy use
- Metal legs extend from 25 cm (10") to 90 cm (36")



ARMS COLLAPSED



ARMS EXTENDED

Reference code	Description	Unit
12-1034HR	XTender <sup>®</sup> extendable goniometer	1 pc



**Lafayette goniometer** MD CE


A goniometer is a useful tool to assess the range of motion of the joints. The Lafayette Goniometer has extendable arms, which allows it to lie over anatomical bony landmarks. This results in the ability to measure all joints, no matter how big or small. The goniometer also has a magnifying lens to read the double 180 degree scale.



Reference code	Description	Unit
12-1035 (old reference: O8-O30701)	Lafayette Gollehan Extendable Goniometer with Magnifying Lens   24 x 26 cm   180°	1 pc

**Plastic goniometers** MD CE


- Scale reads in 1 degree increments; linear scale reads in cm and inches
- Hi-Res<sup>®</sup> option has white background behind the gradations and numbers for enhanced readability
- Take readings from absolute vertical or horizontal with the Absolute+Axis<sup>®</sup> goniometer




CLINIC - 15 CM PLASTIC ISOM (BOTTOM RIGHT)  
CLINIC - 20 CM PLASTIC ISOM (BOTTOM LEFT)  
CLINIC - 30 CM PLASTIC ISOM (TOP)




PATIENT - 15 CM POCKET



30 CM ISOM ABSOLUTE+AXIS<sup>®</sup> (12-1025)



30 CM ISOM ABSOLUTE+AXIS<sup>®</sup> HI-RES<sup>®</sup>  
(INTEGRATED INTO GONIOMETER) (12-1025HR)



ABSOLUTE+AXIS<sup>®</sup> ATTACHMENT SLIDES ONTO  
ANY 30 CM GONIOMETER

Reference code	Description	Unit
Clinic - plastic ISOM		
12-1000HR	30 cm (12") ISOM Hi-Res <sup>®</sup> 360°	1 pc
12-1001HR	20 cm (8") ISOM Hi-Res <sup>®</sup> 360°	1 pc
12-1002HR	15 cm (6") ISOM Hi-Res <sup>®</sup> 360°	1 pc
Patient - plastic pocket		
12-1005HR	15 cm (6") pocket Hi-Res <sup>®</sup> 180°	1 pc
Plastic sets		
12-1028HR	Hi-Res <sup>®</sup> set - 6 pieces	Set
Absolute axis - plastic		
12-1025 (old reference: O8-O30401)	30 cm (12") Absolute+Axis <sup>®</sup>	1 pc
12-1025HR	30 cm (12") Absolute+Axis <sup>®</sup> Hi-Res <sup>®</sup>	1 pc
12-1026	AA attachment only	1 pc




## Digital goniometer

MD CE

- Digital plastic goniometer. Reads 0-360° on a LCD screen; 25 cm linear scale

Reference code	Description	Unit
12-1049	25 cm (10") plastic 360°	1 pc



## Touch-Test® 2-point discriminator

MD CE

- Tests static, **moving** 2-point discrimination
- Consists of two sturdy, rotating, plastic discs that are joined together
- Includes instructions and carrying case

Reference code	Description	Unit
12-1495	Touch Test® 2-point	1 pc




## Discrim-A-Gon® 2-point discriminator


MD CE

- Easy-to-use, lightweight plastic wheel is the perfect sensory evaluation tool to test static and dynamic 1 and 2-point discrimination
- 2 separate 2-point discrimination octagons (D1 and D2) come in a set
- Each octagon measures a different range of 8 labeled, **fixed**, 2-point intervals ranging from 1 to 25mm

Reference code	Description	Unit
12-1492	Discrim-A-Gon® set (2 discs)	Set



D1




D2

## 2-point discriminators

MD CE

- Evaluate cutaneous sensitivity and touch threshold
- Third point permits alternating between single and double point stimulation without changing setting
- Plastic points minimize influence of temperature
- Measures up to 4 cm

Reference code	Description	Unit
12-1480	2-point discriminator	1 pc
12-1481	with third point	1 pc





## Sit 'n Reach® box test

MD CE

- Perform sit 'n reach® flexibility tests
- Measurements include: fitness program, Cooper Institute/YMCA, SHAPE America (formally known as AAHPERD)
- Scales printed in cm and inches
- Sturdy metal construction with a built-in foot plate and a maximum stretch indicator that remains until reset

Reference code	Description	Unit
12-1085	Standard Sit 'n Reach® box test	1 pc



## Finger flexion gauge

MD CE

- Clear plastic gauge simplifies and speeds measurement of composite finger flexion and finger-thumb opposition
- Scale measures up to 10 cm

Reference code	Description	Unit
12-1060	Finger flexion gauge	1 pc




## Gulick measurement tape

MD CE

- Spring exerts constant force on tape for accurate and repeatable results
- Push-button retractor
- Reads in cm and inches

Reference code	Description	Unit
12-1201	150 cm tape	1 pc
12-1203	180 cm tape	1 pc




## Finger circumference gauge

MD CE

- Wrap webbing around finger or digit and read circumference in both cm and inches

Reference code	Description	Unit
12-1222	Gauge (in/cm)	1 pc




## Circumference tape

MD CE

- Measure the circumference of a body part in both cm and inches

Reference code	Description	Unit
12-1205	Tape	1 pc





## Clinical monofilament sensory evaluators



The Baseline Monofilament evaluator is used to measure cutaneous sensory perception threshold. It measures sensitivity, nerve damage or impairment, and changes in neurological status. Each monofilament represents a unique amount of force which is marked by force and size.

### 1. TACTILE™ MONOFILAMENTS

These monofilaments have a fitted protective case that enables safe storage and transport.


Reference code	Description
12-1660 (old reference: O8-O5O421)	2-piece set, normal   2,83 - 3,61
12-1668 (old reference: O8-O5O425)	2-piece set, protective sensation   4,56 - 5,07



### 2. FOLD-UP™ MONOFILAMENTS

These monofilaments feature an easy to fold nylon filament for protection as well as easy storage and portability. Sets are shipped and stored in a soft-shell protective carrying case.

Reference code	Description
12-1740	5-piece set, hand   2,83 - 3,61 - 4,31 - 4,56 - 6,65
12-1741	6-piece set, foot   2,83 - 3,61 - 4,31 - 4,56 - 5,07 - 6,65
12-1743	20-piece set, complete




20-PIECE SET:				
COLOUR	USE	QUANTITY	SIZE	FORCE
●	Normal	1 pc	1,65	(0,008 gm)
●		1 pc	2,36	(0,02 gm)
●		1 pc	2,44	(0,04 gm)
●		1 pc	2,83	(0,07 gm)
●	Diminished light touch	1 pc	3,22	(0,16 gm)
●		1 pc	3,61	(0,4 gm)
●	Diminished protective sensation	1 pc	3,84	(0,6 gm)
●		1 pc	4,08	(1,0 gm)
●		1 pc	4,17	(1,4 gm)
●		1 pc	4,31	(2,0 gm)
●	Loss of protective sensation	1 pc	4,56	(4,0 gm)
●		1 pc	4,74	(6,0 gm)
●		1 pc	4,93	(8,0 gm)
●		1 pc	5,07	(10,0 gm)
●		1 pc	5,18	(15,0 gm)
●		1 pc	5,46	(26,0 gm)
●		1 pc	5,88	(60,0 gm)
●		1 pc	6,10	(100,0 gm)
●		1 pc	6,45	(180,0 gm)
●		1 pc	6,65	(300,0 gm)



## 3. DISPOSABLE MONOFILAMENTS

These are individual 5,07 / 10 gm monofilaments that are disposable, making them ideal for clinics and other multi-patient environments

Reference code	Description	Unit
12-167O	Baseline Tactile Monofilament   LEAP Program   Disposable	1 pc
12-167O-2O	Baseline Tactile Monofilament   LEAP Program   Disposable	2O pcs
12-167O-4O	Baseline Tactile Monofilament   LEAP Program   Disposable	4O pcs



## 4. MULTI-FILAMENT PEN



The multi-filament pen offers a non-invasive evaluation of the diminished protective sensation. The pen contains 4 nylon monofilaments, each representing a unique amount of force which is marked by colour and size. Simply slide down the lever of the desired monofilament and use as intended. This pen provides a useful set of monofilaments in a fun, compact, and protected form.

COLOUR	USE	SIZE	FORCE
●	Normal	2,83	(O,07 gm)
●	Diminished light touch	3,61	(O,4 gm)
●	Diminished protective sensation	4,56	(4 gm)
●	Loss of protective sensation	5,07	(10 gm)

Reference code	Description	Unit
12-137O	Baseline Multi-Filament Pen	1 pc



## Tuning forks



- Use to evaluate audio and vibratory sensation

Reference code	Description	Unit
Set of six pieces with case		
12-146O	128, 256, 512, 1O24, 2O48, 4O96 cps	6 pcs
Weighted tuning forks		
12-1466	128 cps weighted	1 pc
12-1467	256 cps weighted	1 pc
Unweighted tuning forks		
12-1468	512 cps	1 pc




## Pinwheel



- Test for nerve damage or nerve dysfunction with the **Baseline® Clearwheel®**
- Used by orthopedists, neurologists, chiropractors and other healthcare practitioners to test a patient's ability to feel sharp contact with the skin
- Consists of a rotating pinwheel attached to a handle

Reference code	Description	Unit
12-1446	Non-sterile Clearwheel®	1 pc





## Bubble<sup>®</sup> inclinometer MD CE

- Place the inclinometer near the joint to measure, turn dial until scale reads 0, take joint through its range, read range from dial
- Some measurements need the use of 2 inclinometers
- Standards codified in "AMA Guide to the Evaluation of Permanent Impairment, third edition"

Reference code	Description	Unit
12-1056	Bubble <sup>®</sup> inclinometer	1 pc



## AcuAngle<sup>®</sup> inclinometers MD CE

- Place inclinometer near joint to be measured; turn dial to 0; take joint through its range; read ROM from dial
- AcuAngle<sup>®</sup> Inclinometer** has adjustable feet along scale (side-to-side) that adapts to body contours; Pointer is dampened by fluid to assure accurate ROM measurements
- AcuAngle<sup>®</sup> Wrist Inclinometer** measures ROM of the wrist along all axes; distance between feet is adjustable

Reference code	Description	Unit
12-1149	AcuAngle <sup>®</sup> inclinometer	1 pc
12-0502	AcuAngle <sup>®</sup> wrist inclinometer	1 pc



ACUANGLE<sup>®</sup> WRIST INCLINOMETER



ACUANGLE<sup>®</sup> INCLINOMETER

## Digital inclinometer MD CE

- Place inclinometer near the joint to be measured; press "zero" button; move joint through its range; press the "hold" button
- Read ROM directly after the joint has been taken through its range



Reference code	Description	Unit
12-1057	Digital inclinometer	1 pc

## Scoliometer MD CE

- Measure the degree of rotation of a deformity on the back
- Small, non-invasive
- 19 cm (7,5") L x 9 cm (3,5") W
- Includes storage pouch

Reference code	Description	Unit
12-1099	Scoliometer	1 pc





## Economy skinfold caliper

MD CEE

- **Plastic** skinfold caliper is the inexpensive choice for skinfold analysis
- Spring-loaded arm assures reproducible results
- Floating tips assure accurate and reproducible results
- 50 mm scale
- With booklet and tables



WITH FLOATING MEASUREMENT TIPS

Reference code	Description	Unit
12-1112	Economy skinfold caliper	1 pc

## Skinfold calipers

MD CEE

- 2-sided lightweight **aluminum** caliper is easy to use and read for left- or right-handed people; includes calibration block
- Floating tips to assure accurate and reproducible results
- Large, easy-to-read 70 mm scale
- Includes protective carrying case
- Calibration block verifies at 5-50 mm
- Easily fits into your pocket



DOUBLE SIDED CALIPER - FRONT



DOUBLE SIDED CALIPER - BACK



CALIBRATION BLOCK



SINGLE SIDED CALIPER

Reference code	Description	Unit
12-1130	Deluxe (dual sided) with calibration block	1 pc + block
12-1110	Skinfold caliper (single sided) without calibration block	1 pc

## BMI Body Fat Scale

MD CEE

- Made with tempered glass
- Measures weight, body fat and water %
- Measurements range from 4% to 60% body fat; 27.5% to 66% hydration; 20% to 56% muscle; 2% to 20% bone
- Percentages are measured in 0.1% increments
- 12-person memory
- 32 cm (12.5") L x 31 cm (12.25") W x 5 cm (2") H
- Weight capacity is 400 lb / 180 kg
- Requires 2 AAA batteries (included)



Reference code	Description	Unit
12-1190 (old reference: 08-070101)	BMI Body Fat Scale	1 pc



# Private label





# Do you want to work together?

You like our products, but you would like your own brand or your customer's name on one or more of them? **Let's work together!**

In close collaboration with our manufacturers from all the products in our range, we can create your desired product with the exact look you want. We can **replace our logo with yours, change the colour, even the packaging** can be tweaked to your needs and liking. Private labels will allow you to design the product the way you want to match it to your company's or customer's identity.

Private Label products are identical in quality as all our products, with the benefit of putting **your brand in the spotlights**. Looking for the perfect solution to a homogeneous range of products? Want to increase your brand awareness? Need a one time opportunity, such as a promotional product? We'll take care of it.

Please **contact us** for more information about pricing and volume.  
sales@mvs-in-motion.com  
marketing@mvs-in-motion.com



Product description	Minimum order quantity
Resistive Exercise Band	100 units
Resistive Exercise Loop	100 units
Soft-Over-Ball	1000 pcs per colour and dimension
Hollow Foam Roller	300 pcs
AB Gym Ball	5000 pcs per colour and dimension
Comfort Gym Mat	300 pcs per colour
Xtra Comfort Gym Mat	
Pearl Packs	500 pcs
Kinesiology tape	600 pcs per colour
Various other products	Contact us for more info



# After-sales





## MVS In Motion calibration service

Measurement Instruments are calibrated at manufacturing with a full one year guarantee, before sending them to users all over the world. Ideally, these instruments should be recalibrated every year.

In order to make this happen in an easy and inexpensive way, Hydraulic Hand Dynamometers (Saehan, Baseline, ...) and Hydraulic Pinch Gauges are recalibrated by our experts at the MVS In Motion Centre in Roeselare - Belgium. This operation consists of re-establishing the exact relationship between the force applied to a dynamometer and the measurement it shows. Small repairs can be done as well.

Simply send your request to our customer service department for a return authorisation. After each recalibration, you will receive an updated calibration certificate, certifying accurate and repeatable test results.

If you have many calibrations to do, you might consider buying your own calibration unit.



Reference Code	Description	Unit
08-019901	Calibration Unit	1 pc
08-019900	Calibration & Repair	-

## MVS In Motion repair centre

Since machines are designed and assembled by humans, they are not flawless. That is why MVS In Motion installed an inhouse repair centre, to swiftly and correctly solve any problems that may occur with our products.

No matter the cause, be it faulty fabrication or misuse, our technician will shine his light on the problem.

Send an email to our customer service department to inform us on the type of problem and the serial number / production date of the device. After our authorisation, please send over the patient, accompanied by your name and client number and a brief description of the symptoms to [aftersales@mvs-in-motion.com](mailto:aftersales@mvs-in-motion.com). We'll do our best to have it **cured in no time!**





# Terms & Conditions





## Article 1. Company details

These sales conditions apply to all current and future agreements with MVS IN MOTION BV, with registered office at 8800 Roeselare, Krommebeekpark 11 and registered with the RPR Gent, Division Kortrijk under number 0458.496.729 (hereinafter referred to as "MVS"), as well as with any company affiliated to MVS in accordance with Article 1:20 of the Companies and Associations Code.

## Article 2. Applicability

These conditions are applicable to all offers, all agreements, all purchases and orders, both made in MVS's office and online, by e-mail or any other medium (hereinafter generally referred to as "Agreement(s)"). As soon as the customer (hereinafter referred to as "Customer(s)") enters into an Agreement with MVS, it accepts the general conditions as described in this document to the exclusion of the application of any of its own general conditions, even if these contain the same grounds for exclusion. The General Terms and Conditions are applicable to all products and services offered by MVS, except in the case of an explicit deviation which is the result of mutual agreement and written confirmation. Explicit deviations shall only be valid as a replacement or an addition to the clauses to which they relate. The other provisions of these general terms and conditions shall continue to apply in full. If a provision contained in these general terms and conditions is contrary to an international treaty, insofar as it is applicable, a law, regulation or other mandatory regulation, the provision in question shall not apply. The invalidity of a provision shall not affect the validity of the remaining provisions of the general conditions, which shall continue to apply in full. MVS reserves the right to change and/or to supplement the general conditions at any time. An amendment of these general conditions shall have no effect on already existing Agreements. Any matters not explicitly provided for in these general terms and conditions of sale shall be governed by the provisions of Belgian law.

## Article 3. Conclusion of the Agreement

Anyone placing an order confirms that they are commercially competent. Commercially incompetent persons shall be represented by their legal representative, who shall also be bound by these general conditions. MVS undertakes to process orders as long as stocks last and within the limits specified in these general conditions. The order can only be processed when the Customer has clearly identified himself. The placing of an order shall constitute acceptance of these general sales conditions. Placing an order entails an obligation to pay. The general conditions are always available and can be consulted on the website. MVS reserves the right to refuse orders in case of serious suspicion of abuse of rights or bad faith, serious suspicion of commercial purposes unacceptable to MVS or in case of exhaustion of stock of a certain article. For each order, the general terms and conditions that are included on the website and/or on the invoice are applicable. The language of the Agreements concluded is English. An order is only final after acceptance of these general terms and conditions, the prices and the description of the offer. If anything is unclear in the offer or after ordering, the Customer can always contact us, either via [info@mvs-in-motion.com](mailto:info@mvs-in-motion.com) or by telephone between 9 a.m. and 5 p.m. on +32 51 50 55 01.

## Article 4. Signature and proof

The Customer accepts electronic evidence. The Customer is solely responsible for the accuracy of all the data he or she transmits. The final confirmation of the order by the Customer is deemed to be acceptance of the order at the price set.

## Article 5. Prices and Payment

The prices are ex works Roeselare (EXW), Belgium. All prices of products are indicated in Euros and are exclusive of VAT. Shipping costs, transport costs, bank charges, delivery costs, taxes and other levies are not included in MVS' prices and will be paid by the Customer, unless otherwise agreed in writing. For orders with an order value of less than 250 € an administrative cost of 10 € will be added. The Customer owes MVS the price that was indicated in the confirmation of the delivery note. Obvious errors in the price quote, such as evident inaccuracies, can be corrected by MVS even after the conclusion of the agreement. All offers made by MVS are free of obligation. MVS can change prices at any time without prior notification. In the absence of timely or full payment, default interest of twelve percent (12%) per year will be owed, ipso jure and without notice of default, from the due date of the invoice until the day of full payment. In the event of non-payment, late payment or full payment, a fixed compensation amounting to 15% of the invoice amount, costs included, with a minimum of EUR 40 will also be payable. In case of non-timely payment of an invoice, all other claims of the Customer that have not yet fallen due will become immediately due by operation of law. With non-payment, MVS reserves the right not to deliver the products. MVS also reserves the right to regard the sales agreement as dissolved by operation of law and without prior notice of default, for the whole or the part not yet executed.

## Article 6. Delivery of the Order

Products will only be delivered in the countries for which MVS allows delivery. The chosen delivery option influences the delivery period. The stated delivery times are not binding, but are given merely as an indication. Incorrectly communicated delivery addresses are the responsibility of the Customer and may result in additional costs. MVS is entitled to make partial deliveries. The Customer is responsible for all additional taxes and duties. MVS cannot be held responsible for any fees charged by a third party. The delivery will be at the risk of the Customer. No delay in delivery can lead to the termination of the agreement, except in case of deliberate delay. MVS commits to pack the products to be delivered as carefully as possible in order to exclude the possibility of damage during normal handling. MVS accepts no responsibility for late deliveries due to unforeseen circumstances caused by third parties, nor for loss, theft or damage during shipment by third parties and, insofar as permitted by law, MVS shall also not be liable for losses, costs, liability, damages, fines or expenses caused by a delay or failure to deliver due to unforeseen circumstances caused by third parties, nor for loss, theft or damage during shipment by third parties.

## Article 7. Reservation of Ownership

The ownership of products will only be transferred to the Customer after full payment of all sums due by the Customer to MVS, including payment of the agreed price, costs, interest and possible damages. Until the actual transfer of ownership, the Customer is prohibited from reselling, renting out, pledging, processing, incorporating or in any other way disposing of or encumbering the products until all payments and obligations have been fully settled, unless otherwise agreed by MVS. Nevertheless, the Customer shall bear the risk of loss or destruction of the Products from the moment of delivery.

## Article 8. Force Majeure

MVS will not be responsible for delays in the fulfilment or non-fulfilment of its commitments due to events beyond its normal control, including production interruptions, difficulties in acquiring or shortages of raw materials, work force, energy or transport or delays in transport, strikes, epidemics, pandemics, lock-outs, work interruptions or other collective labour disputes, which affect either MVS or its suppliers, disruptions in the internet, the electricity network, mail traffic or technology provided by third parties, government measures and illness, even if these events are foreseeable.

## Article 9. Intellectual Property Rights

All intellectual property rights and derived rights are retained by MVS. These intellectual property rights include copyrights, trademark rights, drawing and model rights and/or other (intellectual property) rights, including patentable or non-patentable technical and/or commercial know-how, methods and concepts. The contents of the website, images, logos, drawings, photos, data, product names, texts, etc. are protected by copyright and other (intellectual) property rights and treaty provisions, which exclusively belong to MVS. The photographic work of MVS is subject to copyrights, which are exclusively retained by MVS. For use for commercial purposes, the Customer must request permission and/or pay a fee to MVS.

## Article 10. Processing of Customer Data and Privacy

When products are purchased, the Customer explicitly agrees that his or her personal data can be used and processed for administrative purposes, such as keeping track of orders in a customer database or for billing purposes. The information provided by the Customer is necessary for the processing and completion of orders and invoicing. If these data are missing, the order will inevitably be cancelled. Personal data will only be processed for promotional purposes, personal advertising and/or other (re)marketing purposes if the Customer has given his/her consent. MVS will never pass on this personal information to third parties, unless with the Customer's consent or when explicitly requested by judicial authorities or when it is legally obliged to do so. The Customer has the right to modify or review his or her personal information at any time. The Customer also has the right to oppose, free of charge, the processing of his personal data for direct marketing initiatives. Only Customers who have given their express consent to this will be included in MVS's newsletter. Unsubscribing is possible upon simple request of the Customer. Should the Customer wish to contact MVS with regard to the processing of his or her personal data, a message should be sent to the following e-mail address: [info@mvs-in-motion.com](mailto:info@mvs-in-motion.com)

## Article 11. Defects and complaints

The Customer undertakes to examine the products thoroughly immediately upon receipt. MVS' warranty is limited to the warranty for hidden defects and the warranty in the case of consumer purchases, both as specified in the Civil Code. MVS is only responsible for damages which are the result of intent or gross error. MVS is not responsible for general or special indirect damages of any kind suffered by the Customer nor for damages to goods, visible at delivery, which are not mentioned on the carrier's consignment note. MVS shall also not be responsible for any postponement or cancellation of a production that is due to late delivery by the Customer or its subcontractor of a product or accessory, or to the punctuality or faulty conduct of the latter. Complaints must reach MVS in order to be admissible: (a) in the case of complaints due to non-conformity of the delivery, seven (7) days after delivery and (b) in the case of hidden defects, seven (7) days after discovery of the defect or after the moment by which such discovery could reasonably have taken place. The acceptance of the goods by the Customer covers any possible defect or non-conformity which could have been discovered at that time. MVS undertakes to answer the complaint within a period of 14 calendar days after written receipt. Reimbursement or exchange are only possible after explicit agreement by MVS. Prior to the return of the products, the Customer will provide MVS with a detailed inventory of the products that will be returned (cfr Complaint Form) as well as any other information that MVS may require, including but not limited to the batch number, the production date or the expiry date. Under no circumstances will returns of more than 30 days after delivery be accepted, nor will returns be accepted for products that are damaged, used or not in their original and unopened packaging. In the event of a return, MVS may charge the Customer a processing fee of ten percent (10%) of the purchase price (costs and taxes included).

## Article 12. Applicable law and competent court

All offers and Agreements are governed exclusively by Belgian law. In the event of disputes, the courts of the judicial district of Gent shall have jurisdiction. Within the framework of their relations, the parties concerned accept electronic evidence.









**MVS In Motion** | Krommebeekpark 11, 8800 Roeselare | Belgium - Europe  
+32 51 50 55 01 | [www.mvs-in-motion.com](http://www.mvs-in-motion.com) | [info@mvs-in-motion.com](mailto:info@mvs-in-motion.com)